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
Created on August 1, 2020

## Background

Angiotensin-converting enzyme 2 (ACE2), also called ACE homolog, is a membrane protein considerably homologous to ACE.

ACE2 was known as a factor in the renin-angiotensin system involved in blood pressure regulation.

Recently, it was clarified that the Spike glycoprotein of SARS-CoV-2 binds to ACE2 present on the cell membrane and then the virus is taken up into the cell when the infection occurs. Therefore, ACE2 is considered as a receptor for SARS-CoV-2 <sup>(1)</sup>.


As shown in the usage example, a highly sensitive binding assay can be performed by this product (HAK-ACE2\_UL-1) and the receptor binding domain of SARS-CoV-2 Spike glycoprotein (HAK-SPD\_UL-1, HAK-SPD\_bio-1).

The product must be stored at -70°C or lower upon receipt.

It is recommended to make small aliquots when using.

## 1. SDS-PAGE

430ng of the product was electrophoresed on 4~20% gradient gel and stained with CBB.

## Product information

### 1. Protein Construction

Human ACE2, His Tag (HAK-ACE2\_UL-1) contains Gln18-Pro738 of human ACE2 protein tagged with polyhistidine residues at its C-terminus. The protein was expressed and secreted from HEK293 cells, followed by the purification using a Ni column.

### 2. Purity

> 95% (SDS-PAGE)

### 3. Composition

Human ACE2, His Tag 1mg/mL,  
0.1M-PBS(pH7.2~7.4)


### 4. Storage:

Please avoid repeated freeze-thaw cycles.

## Usage example

### 1. Binding assay

Human ACE2 His tag was immobilized in a 96-well plate by adding at serial dilution of 0.156-10µg/mL. Then 0.8µg/mL of biotin labeled SARS-CoV-2


Spike glycoprotein Receptor-binding

domain, Rabbit IgG1 Fc Tag was added to the wells. The binding was detected with HRP-labeled streptavidin. The linear range in the standard curve was 0.156-1.25µg/mL. Detail protocol for the binding assay is provided separately.

## 2. Screening for the binding inhibitors by competition assay

Human ACE2 His tag was immobilized in a 96-well plate by adding at 1µg/mL. Then 600ng/mL biotin labeled SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc tag and a binding inhibitor (SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc tag in this case) were added to the well and the binding was detected by HRP-labeled streptavidin. The IC<sub>50</sub> was 140ng/mL (2.7nM).

## References

- (1) Alexandra C. Walls, Young-Jun Park, et al., Cell 180, 1-12 (2020)

# Protocol for the binding assay of SARS-CoV-2 Spike glycoprotein and ACE2

Code No. HAK-SPD\_bio-1

Code No. HAK-SPD\_UL-1

Code No. HAK-ACE2\_UL-1

August 1, 2020

---

## 【Experiment 1, binding assay】

### Reagents


- Human ACE2, His Tag
- Biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
- HRP-labeled streptavidin
- Immobilization buffer : PBS (-)
- Dilution buffer/Blocking buffer : 1%BSA/ PBS (-)
- Washing buffer : PBS(-) containing 0.05%Tween20
- HRP substrate : TMB etc.
- Stopping solution : 2N H<sub>2</sub>SO<sub>4</sub>
  
- 96-well plate
- Plate seal
- Plate shaker
- Plate reader

### Preparation of reagents

1. Dilute human ACE2, His Tag with PBS(-) to 1 μg /mL. 100 μL of the solution is needed per well.
  
2. Dilute biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 29-21000ng/mL (serial dilution by three-fold). 100 μL of the solution is needed per well.
  
3. Dilute HRP-labeled streptavidin according to manufacturer's manual. 100 μL of the solution is needed per well.

### Method for the measurement

1. Dilute human ACE2, His Tag.
2. Add 100  $\mu\text{L}$ /well of human ACE2 in 96-well plate.
3. Stand at 4°C overnight.
4. Remove the reagent completely and add 300 $\mu\text{L}$ /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
5. Add 200 $\mu\text{L}$  /well of 1%BSA/PBS(-).
6. As a blocking procedure, stand the plate at room temperature for 2 hours or at 4°C overnight.
7. Remove the reagent completely and add 300 $\mu\text{L}$ /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
8. Add 100 $\mu\text{L}$ /well of diluted biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag.
9. Stand at room temperature for 2 hours.
10. Remove the reagent completely and add 300 $\mu\text{L}$ /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
11. Add 100 $\mu\text{L}$ /well of diluted HRP-labeled streptavidin.
12. Stand at room temperature (reaction time should be referred to the manufacturer's manual) .
13. Remove the reagent completely and add 300 $\mu\text{L}$ /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
14. Add 100 $\mu\text{L}$ /well of substrate solution to react at room temperature.
15. When the color developed, add 50 $\mu\text{L}$ /well of stopping solution.
16. Measure the optical density (OD450) of each well by plate reader.

## 【Experiment 2 : Screening for the binding inhibitors by competition assay】

### Reagents

- Human ACE2, His Tag
  - Bioin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
  - SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
  - HRP-labeled streptavidin
  - Immobilization buffer : PBS (-)
  - Dilution buffer/Blocking buffer : 1%BSA/ PBS (-)
  - Washing buffer : PBS(-) containing 0.05%Tween20
  - HRP substrate : TMB etc.
  - Stopping solution : 2N H<sub>2</sub>SO<sub>4</sub>
- 
- 96-well plate
  - Plate seal
  - Plate shaker
  - Plate reader

### Preparation of reagents

1. Dilute human ACE2, His Tag with PBS(-) to 1μg /mL. 100 μL of the solution is needed per well.
2. Dilute biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 1200ng/mL. 50 μL of the solution is needed per well to make final concentration of 600ng/mL.
3. Dilute SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 20-10000ng/mL (serial dilution by two-fold) to make final concentration of 10-5000ng/mL.
4. Dilute HRP-labeled streptavidin according to manufacturer's manual.

### Method for the measurement

1. Dilute human ACE2, His Tag.
2. Add 100 μL/well of human ACE2 in 96-well plate.
3. Stand at 4°C overnight.
4. Remove the reagent completely and add 300μL/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
5. Add 200μL /well of 1%BSA/PBS(-).
6. As a blocking procedure, stand the plate at room temperature for 2 hours or at 4°C

overnight.

7. Remove the reagent completely and add 300 $\mu$ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
8. Add 50 $\mu$ L/well of diluted biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag.
9. Add 50 $\mu$ L/well of SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag as a control inhibitor or any inhibitor samples to be tested.
10. Stand at room temperature for 2 hours.
11. Remove the reagent completely and add 300 $\mu$ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
12. Add 100 $\mu$ L/well of diluted HRP-labeled streptavidin.
13. Stand at room temperature (reaction time should be referred to the manufacturer's manual).
14. Remove the reagent completely and add 300 $\mu$ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
15. Add 100 $\mu$ L/well of substrate solution to react at room temperature.
16. When the color developed, add 50 $\mu$ L/well of stopping solution.
17. Measure the optical density (OD450) of each well by plate reader.
18. Check for the decrease in the absorbance caused by the binding inhibitors.