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
Created on August 1, 2020

Background

Angiotensin-converting enzyme 2 (ACE2), also called ACE homolog, is a membrane protein considerably homologous to ACE.

ACE2 was known as a factor in the renin-angiotensin system involved in blood pressure regulation.

Recently, it was clarified that the Spike glycoprotein of SARS-CoV-2 binds to ACE2 present on the cell membrane and then the virus is taken up into the cell when the infection occurs. Therefore, ACE2 is considered as a receptor for SARS-CoV-2 and this product (SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag) corresponds to the binding site of the Spike glycoprotein to the viral receptor ACE2 ⁽¹⁾.

As shown in the usage example, a highly sensitive binding assay can be performed by this product (HAK-SPD_bio-1) and human ACE2 His tag (HAK-ACE2_UL-1).

Product information

1. Protein Structure

Biotin-labeled SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag (HAK-SPD_bio-1) contains Arg319-Phe541 of SARS-CoV-2 Spike glycoprotein tagged with Fc region of rabbit IgG1 at its C-terminus. The protein was expressed and secreted from HEK293 cells, followed by the purification using a Protein A column and then biotinylated.

2. Purity

> 95% (SDS-PAGE)

3. Composition

SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag 1mg/mL

0.1M PBS(pH7.2~7.4)

0.1% BSA

4. Storage:


Please avoid repeated freeze-thaw cycles.

The product must be stored at -70°C or lower upon receipt.

It is recommended to make small aliquots when using.

5. SDS-PAGE

500ng of the product before biotinylation was electrophoresed on 4~20% gradient gel and stained with CBB.

Usage example

1. Binding assay

For research use only, not for diagnostic use.
Please read this manual carefully before use.


Biotin labeled SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag (20µg)

Code No. HAK-SPD_bio-1

Created on August 1, 2020

Human ACE2 His tag was immobilized in a 96-well plate by adding at 0.1µg/100 µL/well. Then 29-21000ng/mL of biotin labeled SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag was added to the wells. The binding was detected with HRP-labeled streptavidin. The linear range in the standard curve was 29-798ng/mL. Detail protocol for the binding assay is provided separately.

For research use only, not for diagnostic use.
Please read this manual carefully before use.

Protocol for the binding assay of SARS-CoV-2 Spike glycoprotein and ACE2

Code No. HAK-SPD_bio-1

Code No. HAK-SPD_UL-1

Code No. HAK-ACE2_UL-1

August 1, 2020

【Experiment 1, binding assay】

Reagents

- Human ACE2, His Tag
- Biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
- HRP-labeled streptavidin
- Immobilization buffer : PBS (-)
- Dilution buffer/Blocking buffer : 1%BSA/ PBS (-)
- Washing buffer : PBS(-) containing 0.05%Tween20
- HRP substrate : TMB etc.
- Stopping solution : 2N H₂SO₄

- 96-well plate
- Plate seal
- Plate shaker
- Plate reader

Preparation of reagents


1. Dilute human ACE2, His Tag with PBS(-) to 1 μg /mL. 100 μL of the solution is needed per well.

2. Dilute biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 29-21000ng/mL (serial dilution by three-fold). 100 μL of the solution is needed per well.

3. Dilute HRP-labeled streptavidin according to manufacturer's manual. 100 μL of the solution is needed per well.

Method for the measurement

1. Dilute human ACE2, His Tag.
2. Add 100 μL /well of human ACE2 in 96-well plate.
3. Stand at 4°C overnight.
4. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
5. Add 200 μL /well of 1%BSA/PBS(-).
6. As a blocking procedure, stand the plate at room temperature for 2 hours or at 4°C overnight.
7. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
8. Add 100 μL /well of diluted biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag.
9. Stand at room temperature for 2 hours.
10. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
11. Add 100 μL /well of diluted HRP-labeled streptavidin.
12. Stand at room temperature (reaction time should be referred to the manufacturer's manual) .
13. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
14. Add 100 μL /well of substrate solution to react at room temperature.
15. When the color developed, add 50 μL /well of stopping solution.
16. Measure the optical density (OD450) of each well by plate reader.

【Experiment 2 : Screening for the binding inhibitors by competition assay】

Reagents

- Human ACE2, His Tag
 - Bioin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
 - SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
 - HRP-labeled streptavidin
 - Immobilization buffer : PBS (-)
 - Dilution buffer/Blocking buffer : 1%BSA/ PBS (-)
 - Washing buffer : PBS(-) containing 0.05%Tween20
 - HRP substrate : TMB etc.
 - Stopping solution : 2N H₂SO₄
-
- 96-well plate
 - Plate seal
 - Plate shaker
 - Plate reader

Preparation of reagents

1. Dilute human ACE2, His Tag with PBS(-) to 1μg /mL. 100 μL of the solution is needed per well.
2. Dilute biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 1200ng/mL. 50 μL of the solution is needed per well to make final concentration of 600ng/mL.
3. Dilute SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 20-10000ng/mL (serial dilution by two-fold) to make final concentration of 10-5000ng/mL.
4. Dilute HRP-labeled streptavidin according to manufacturer's manual.

Method for the measurement

1. Dilute human ACE2, His Tag.
2. Add 100 μL/well of human ACE2 in 96-well plate.
3. Stand at 4°C overnight.
4. Remove the reagent completely and add 300μL/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
5. Add 200μL /well of 1%BSA/PBS(-).
6. As a blocking procedure, stand the plate at room temperature for 2 hours or at 4°C

overnight.

7. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
8. Add 50 μ L/well of diluted biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag.
9. Add 50 μ L/well of SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag as a control inhibitor or any inhibitor samples to be tested.
10. Stand at room temperature for 2 hours.
11. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
12. Add 100 μ L/well of diluted HRP-labeled streptavidin.
13. Stand at room temperature (reaction time should be referred to the manufacturer's manual).
14. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
15. Add 100 μ L/well of substrate solution to react at room temperature.
16. When the color developed, add 50 μ L/well of stopping solution.
17. Measure the optical density (OD450) of each well by plate reader.
18. Check for the decrease in the absorbance caused by the binding inhibitors.