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
Created on August 1, 2020

Background

Angiotensin-converting enzyme 2 (ACE2), also called ACE homolog, is a membrane protein considerably homologous to ACE.

ACE2 was known as a factor in the renin-angiotensin system involved in blood pressure regulation.

Recently, it was clarified that the Spike glycoprotein of SARS-CoV-2 binds to ACE2 present on the cell membrane and then the virus is taken up into the cell when the infection occurs. Therefore, ACE2 is considered as a receptor for SARS-CoV-2 and this product (SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag) corresponds to the binding site of the Spike glycoprotein to the viral receptor ACE2⁽¹⁾.

As shown in the usage example, a highly sensitive binding assay can be performed by this product (HAK-SPD_UL-1) and human ACE2 His tag (HAK-ACE2_UL-1).

Product information

1. Protein Structure

SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag (HAK-SPD_UL-1) contains Arg319-Phe541 of SARS-CoV-2 Spike glycoprotein tagged with Fc region of rabbit IgG1 at its C-terminus. The protein was expressed and secreted from HEK293 cells, followed by the purification using a Protein A column.

2. Purity

> 95% (SDS-PAGE)

3. Composition

SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag 1mg/mL,
0.1M PBS(pH7.2~7.4)

4. Storage:


Please avoid repeated freeze-thaw cycles.

The product must be stored at -70°C or lower upon receipt.

It is recommended to make small aliquots when using.

5. SDS-PAGE

500ng of the product was electrophoresed on 4~20% gradient gel and stained with CBB.

Usage example

For research use only, not for diagnostic use.

Please read this manual carefully before use.

SARS-CoV-2 Spike glycoprotein Receptor-

binding domain, Rabbit IgG1 Fc Tag (100µg)


Code No. HAK-SPD_UL-1

Created on August 1, 2020

1. Binding assay

Human ACE2 His tag was immobilized in a 96-well plate by adding at 0.1µg/100 µL/well. Then 0.293-1200ng/mL of SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag was added to the wells. The binding was detected with HRP-labeled anti-rabbit IgG secondary antibody. The linear range in the standard curve was 0.293-18.8ng/mL. Detail protocol for the binding assay is provided separately.

References

- (1) Alexandra C. Walls, Young-Jun Park, et al., Cell 180, 1-12 (2020)

For research use only, not for diagnostic use.
Please read this manual carefully before use.

Protocol for the binding assay of SARS-CoV-2 Spike glycoprotein and ACE2

Code No. HAK-SPD_bio-1

Code No. HAK-SPD_UL-1

Code No. HAK-ACE2_UL-1

August 1, 2020

【Experiment 1, binding assay】

Reagents

- Human ACE2, His Tag
- Biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
- HRP-labeled streptavidin
- Immobilization buffer : PBS (-)
- Dilution buffer/Blocking buffer : 1%BSA/ PBS (-)
- Washing buffer : PBS(-) containing 0.05%Tween20
- HRP substrate : TMB etc.
- Stopping solution : 2N H₂SO₄


- 96-well plate
- Plate seal
- Plate shaker
- Plate reader

Preparation of reagents

1. Dilute human ACE2, His Tag with PBS(-) to 1 μg /mL. 100 μL of the solution is needed per well.
2. Dilute biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 29-21000ng/mL (serial dilution by three-fold). 100 μL of the solution is needed per well.
3. Dilute HRP-labeled streptavidin according to manufacturer's manual. 100 μL of the solution is needed per well.

Method for the measurement

1. Dilute human ACE2, His Tag.
2. Add 100 μL /well of human ACE2 in 96-well plate.
3. Stand at 4°C overnight.
4. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
5. Add 200 μL /well of 1%BSA/PBS(-).
6. As a blocking procedure, stand the plate at room temperature for 2 hours or at 4°C overnight.
7. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
8. Add 100 μL /well of diluted biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag.
9. Stand at room temperature for 2 hours.
10. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
11. Add 100 μL /well of diluted HRP-labeled streptavidin.
12. Stand at room temperature (reaction time should be referred to the manufacturer's manual) .
13. Remove the reagent completely and add 300 μL /well of washing buffer. Discard the washing buffer and repeat this washing process twice.
14. Add 100 μL /well of substrate solution to react at room temperature.
15. When the color developed, add 50 μL /well of stopping solution.
16. Measure the optical density (OD450) of each well by plate reader.

【Experiment 2 : Screening for the binding inhibitors by competition assay】

Reagents

- Human ACE2, His Tag
 - Bioin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
 - SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag
 - HRP-labeled streptavidin
 - Immobilization buffer : PBS (-)
 - Dilution buffer/Blocking buffer : 1%BSA/ PBS (-)
 - Washing buffer : PBS(-) containing 0.05%Tween20
 - HRP substrate : TMB etc.
 - Stopping solution : 2N H₂SO₄
-
- 96-well plate
 - Plate seal
 - Plate shaker
 - Plate reader

Preparation of reagents

1. Dilute human ACE2, His Tag with PBS(-) to 1μg /mL. 100 μL of the solution is needed per well.
2. Dilute biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 1200ng/mL. 50 μL of the solution is needed per well to make final concentration of 600ng/mL.
3. Dilute SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag with 1%BSA/PBS(-) to 20-10000ng/mL (serial dilution by two-fold) to make final concentration of 10-5000ng/mL.
4. Dilute HRP-labeled streptavidin according to manufacturer's manual.

Method for the measurement

1. Dilute human ACE2, His Tag.
2. Add 100 μL/well of human ACE2 in 96-well plate.
3. Stand at 4°C overnight.
4. Remove the reagent completely and add 300μL/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
5. Add 200μL /well of 1%BSA/PBS(-).
6. As a blocking procedure, stand the plate at room temperature for 2 hours or at 4°C

overnight.

7. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
8. Add 50 μ L/well of diluted biotin labeled SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag.
9. Add 50 μ L/well of SARS-CoV-2 Spike glycoprotein, Rabbit IgG1 Fc tag as a control inhibitor or any inhibitor samples to be tested.
10. Stand at room temperature for 2 hours.
11. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
12. Add 100 μ L/well of diluted HRP-labeled streptavidin.
13. Stand at room temperature (reaction time should be referred to the manufacturer's manual).
14. Remove the reagent completely and add 300 μ L/well of washing buffer. Discard the washing buffer and repeat this washing process twice.
15. Add 100 μ L/well of substrate solution to react at room temperature.
16. When the color developed, add 50 μ L/well of stopping solution.
17. Measure the optical density (OD450) of each well by plate reader.
18. Check for the decrease in the absorbance caused by the binding inhibitors.