## **IST Fluolid-W Protein/Antibody Labeling Kits**

**IST Fluolid-W Orange 600 Protein/Antibody Labeling Kit**

**IST Fluolid-W Yellow 540 Protein/Antibody Labeling Kit**

**IST Fluolid-W Green 520 Protein/Antibody Labeling Kit**

### **Contents**

IST Fluolid-W succinimidyl ester in DMSO	240 $\mu$ l	× 3
0.2 M Sodium bicarbonate buffer (pH8.3)	480 $\mu$ l	× 3

### **Application**

Labeling of protein or antibody (1.0 mg)

### **Characteristics**

- 1. high quantum yield in the solid state**
- 2. no photobleach**
- 3. high stability for light, heat and pH**
- 4. labeling rate is higher than that of traditional dye**

### **Attention**

Product is research use only.

The fluorescence unit of "IST Fluolid" has high stability under various conditions.


To keep the activity of the succinimidyl ester group, please store at 2-8 °C.

Contents are subject to change without notice.

### **Remarks**

The kits have been optimized for labeling 1.0 mg of an IgG antibody.

Other proteins can also be used at the same protocol. However, in the case of under- or over- labeling, it is necessary to investigate labeling condition.

### **Manufacture**

International Science Technology CO., LTD.  
Kasuga Lab.

### **Distributor**

COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618

# IST Fluolid-W Protein/Antibody Labeling Kits

IST Fluolid-W Orange 600 Protein/Antibody Labeling Kit, cat# ISU-IST004

IST Fluolid-W Yellow 540 Protein/Antibody Labeling Kit, cat# ISU-IST005

IST Fluolid-W Green 520 Protein/Antibody Labeling Kit, cat# ISU-IST006

## Contents

IST Fluolid-W succinimidyl ester in DMSO	240 $\mu$ l
	× 3
0.2 M Sodium bicarbonate buffer (pH8.3)	480 $\mu$ l
	× 3

## Attention

Product is research use only.

The fluorescence unit of "IST Fluolid" has high stability under various conditions.

To keep the activity of the succinimidyl ester group, please store at 2-8 °C.

Contents are subject to change without notice.

## Protocol

The kits have been optimized for labeling 1.0 mg of an IgG antibody.

Other proteins can also be used at the same protocol. However, in the case of under- or over- labeling, it is necessary to investigate labeling condition.

It may be usable under the low DMSO concentration. In the case of the precipitations, they typically have no effect on the reaction. However, adjustments of the protocol may be necessary for yield loss.

Take care of yield loss, when the conditions of the protocol are changed. In this instance, please investigate the conditions.

### (1) Preparation of Protein Solution

Dissolve the 1.0 mg protein in 400  $\mu$ l sodium bicarbonate buffer (0.2 M, pH8.3).

### (2) Reaction Mixture (Labeling Reaction)

Mix the 200  $\mu$ l IST Fluolid-W succinimidyl ester/DMSO solution with the protein solution prepared in Section (1). The reaction mixture is stirred slowly for 1-2 hours at room temperature. In the case of the precipitations, they have no effect on the reaction.

### (3) Purification of Labeled Protein

The labeled protein can be purified from the reaction mixture by size-exclusion purification.

The following purification method can be used.

- GE healthcare/amershambioscience Sephadex series (G-25, G-50, etc.)
- Bio-Rad Bio Gel P Series (P-6, P-30, etc.)
- Millipore Amicon series etc.

If there is an unknown point, please contact us in the following.


Figure 1. Protocol for Protein Labeling with Protein Labeling kits.

## Characteristics

Fluorescent Dye	$\lambda_{abs}/nm^*$	$\lambda_{em}/nm^*$	$\epsilon/M^{-1}cm^{-1}^+$	Molecular Weight
Fluolid-W Orange 600	440	602	15,000	570
Fluolid-W Yellow 540	410	541	12,000	540
Fluolid-W Green 520	395	522	10,000	510

\* Absorbance and fluorescence emission maxima.

+ Extinction coefficient at  $\lambda_{abs}$ .

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