



## Acidic Mucopolysaccharide Assay Kit

### Introduction

Chondrocytes produce extracellular matrix such as collagen, mucopolysaccharide and so on. There are conventional methods to measure acidic mucopolysaccharides using radioisotopes but the area and the amount used regarding radioisotopes are restricted. On the other hand, HPLC analysis is complicated to maintain the system and to extra for making samples. Acidic Mucopolysaccharide Assay Kit (PMC-AK03-COS) provides a convenient system for mucopolysaccharides visualization using Stains All, which combines with acidic mucopolysaccharides, with ease. Stains All normally combines with acidic substance, however the Stains All in this kit stains only acidic mucopolysaccharides of chondrocytes selectively.

### Components/Storage

Components	Quantity	Stability
Staining Solution	4 mL	4-10°C
Buffer	130 mL	4-10°C
Enzyme Reagent	5 vials	4-10°C
Standard Stock Solution (100 ug/mL Chondroitin Sulfate Calibrator)	2 mL	4-10°C

**One kit contains regents for 100 (20 samples × 5 times) samples**

**This kit cannot quantify various glycosaminoglycans individually.**

**This kit cannot quantify the mucopolysaccharides in culture medium.**

### Additional Materials Required

- distilled water

### Prepare reagents

- Reaction Buffer

Warm Staining Solution and Buffer to room temperature. Mix 0.4 mL of Staining Solution and 12.6 mL of Buffer.

**Note:** After mixing, the solution is blue but the color will fade in several minutes. This working solution cannot be stored.

- Standard Solution

Dilute the Standard Stock Solution in distilled water to 100, 50, 25, 12.5 and 0 ug/mL.

**Note:** This working solution can be stored frozen at -20°C.



COSMO BIO Co., LTD.

Inspiration for Life Science

• Enzyme Solution

Add 10 mL of distilled water to one vial of Enzyme Reagent and mix. Filter the solution with a membrane filter.

**Note:** This working solution cannot be stored.

## Protocol

1. Prepare cultured chondrocyte layers or spheroid cell clusters.

**Note:** Using cell clusters is recommended. Refer to culture spheroid cell clusters method to know how to prepare spheroid cell clusters.

2. Remove culture medium and add 0.5 mL of Enzyme Solution. Digest the cells by incubating at 60°C for 1 hour.


3. Add 100  $\mu$ L of digested cell samples or 100  $\mu$ L of Standard Solution to tube, and add 1.3 mL of Reaction Buffer to each tube, and then leave the tubes.

4. The solution will become blue in a few minutes from the addition. Measure the absorption at 650 nm after the addition of 10 to 20 minutes.

**Note:** When the concentration of acidic mucopolysaccharides is more than 120  $\mu$ g/mL, acidic mucopolysaccharides will precipitate. In the case, the sample must be diluted to obtain a concentration of less than 100  $\mu$ g/mL. After the addition more than 20 minutes, acidic mucopolysaccharides will form blue precipitation

**Fig.1 Acidic Mucopolysaccharide and acidic substances detected by Acidic Mucopolysaccharide Assay Kit**

This kit stains only acidic mucopolysaccharides of chondrocytes selectively.



*For research use only. Not for clinical diagnosis.*



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617


[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619



## Culture Spheroid Cell Clusters Method



**Fig.2 Cells form spheroid clusters under Hypoxia condition**

Prepare cultured chondrocytes and cells are detached by the treatment of trypsin solution. Transfer detached cells to 15 mL centrifuge tube and centrifuge the centrifuge tube. Remove supplement and add 5 mL of medium to the centrifuge tube. In about 3 weeks, cells will form cartilage-like spheroids.

## References

- (1) Obayashi, K., Miyagawa-Tomita, S., Matsumoto, H., Koyama, H., Nakanishi, T., Hirose, H. Effects of Transforming Growth Factor-beta3 and Matrix Metalloproteinase-3 on The Pathogenesis of Chronic Mitral Valvular Disease in Dogs. Am. J. Vet. Res. 72, 194-202 (2011)

## Notice to Purchaser

This product is to be used for Research Purposes Only. It is not to be used for Drug or Diagnostic Purposes, nor is it intended for Human Use. Cosmo Bio products may not be resold, modified for resale, or used to manufacture commercial products without the express written consent of Cosmo Bio Co., Ltd.

Except as otherwise expressly set forth in this use manual, Cosmo Biol dose not make any representation or warranties or conditions of any kind, either express or implied, with respect to the products, or information disclosed hereunder, including, but not limited to, the implied warranties of merchantability, fit for a particular purpose, or noninfringement of the intellectual property rights of third parties.

*For research use only. Not for clinical diagnosis.*



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619