For Research Use Only. Not for use in diagnostic procedures.

c-Myc-tagged Protein

Magnetic PURIFICATION KIT

(mAb. clone 1G4)

CODE No. 3340

*PURIFICATION to Maintain Protein Activity
from eukaryote cell lysate and culture supernatant*

Product Description

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more.

The method of purification is one of the important keys for maintaining protein structure and function. The c-Myc-tagged Protein Magnetic PURIFICATION KIT is designed for the isolation of c-Myc tagged protein from cell lysates and culture supernatant. Severe conditions such as acidic or alkaline elution can denature protein structure while maintaining a neutral pH can preserve protein activity and conformation. MBL has developed the Anti-c-Myc-tag Magnetic Beads to quickly and efficiently purify c-Myc-tagged proteins at neutral pH to maintain protein activity and native conformation. The c-Myc-tag peptide competitively elutes the c-Myc tagged protein from Anti-c-Myc-tag Magnetic Beads while minimally interfering with protein function.

The affinity of the anti-c-Myc-tag antibody is also very important. Very high affinity antibody cannot be dissociated from the c-Myc-tagged protein, and low affinity antibody is not sufficient for binding. MBL has optimized the affinity of the anti-c-Myc-tag antibody to permit high yields and efficient purification of functional c-Myc-tagged proteins for biochemical studies. This kit is simple and fast, because all procedures of this kit have been optimized.

Kit Components.

1. Anti-c-Myc-tag Magnetic Beads 1 mL x 2 vials

10 mg magnetic beads in 1 mL PBS/0.1% BSA with 0.09% sodium azide as preservative

*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.


2. Elution Peptide 1 mg x 2 vials (lyophilized)

c-Myc-tag peptide (EQKLISEEDL), 1 mg in 1 mL PBS after reconstitution

3. Wash Concentrate 5 mL x 2 vials

20 x concentrate

Procedure Summary

Product Capacity

The purification capacity of the Anti-c-Myc-tag Magnetic Beads varies depending upon the c-Myc-tagged protein. For examples, 100 μ L of Anti-c-Myc-tag Magnetic Beads (100 μ L slurry) can be purified 2.2 μ g of a c-Myc-tagged protein (27 kDa).

Materials Required but not Provided

1. Microcentrifuge capable of 15,000 x g
2. Sampling tube (1.5 mL)
3. End-over-end rotator
4. Magnetic Rack (MBL; code 3190)
5. Lysis buffer

Suitable Lysis buffer varies with cell type.

Note: see Additional Information

Homemade Lysis buffer

20-50 mM	Tris-HCl (pH 7.5) or HEPES-KOH (pH 7.5)
50-250 mM	NaCl
5 mM	EDTA
1%	NP-40 or Triton X-100
	if necessary add Protease Inhibitor Cocktail (e.g. SIGMA; code P8340, PIERCE; code 78415).

Storage

Store for up to 1 year from date of receipt at 2-8°C. Do not freeze.

The descriptions of the following protocols are examples. Each user should determine the appropriate condition.

Protocols

The following protocol is for the isolation of c-Myc-tagged proteins produced in a 100-mm cell culture dish. The expression level of the c-Myc-tagged protein may vary. If necessary, adjust the volume of Anti-c-Myc-tag Magnetic Beads and Elution Peptide Solution proportionally.

Material Preparation

1. Wash Solution

Dilute Wash Concentrate with 20 times its volume of distilled water. (e.g. Dilute 0.5 mL of Wash Concentrate with 9.5 mL of distilled water.) For each purification, prepare 5 mL of Wash Solution.

2. Elution Peptide Solution

Reconstitute the Elution Peptide with 1 mL of distilled water. If you want to store the reconstituted Elution Peptide, prepare appropriate aliquots (e.g. 45 μ L x 20 tubes) and store at -20°C. Repeated freezing and thawing is not recommended.

A. Purification from mammalian cultured cell lysate

(c-Myc-tagged protein is not secreted from the cells)

(Lysis of Mammalian Cells)

1. Detach the cells from the culture dish if necessary, and collect the cell suspension into the centrifuge tube.
2. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
3. Wash cells by resuspending the cell pellet in ice-cold PBS.
4. Centrifuge the cell suspension at 400 x g for 5 minutes to pellet the cells. Carefully remove and discard the supernatant.
5. Add 0.5 mL of Lysis buffer to the cell pellet and vortex.
6. Incubate the sample for 15 minutes on ice.
7. Remove cell debris by centrifugation at 15,000 x g for 5 minutes at 4°C.

(Purification of c-Myc-tagged Protein)

8. Transfer the 0.5 mL of cell lysate to the 1.5 mL sampling tube.
9. Resuspend Anti-c-Myc-tag Magnetic Beads by vortexing for a few seconds immediately before dispensing.
10. Dispense 50 μ L of the Anti-c-Myc-tag Magnetic Beads suspension into the sampling tube prepared in step 8. If necessary, please use 100 μ L of the Anti-c-Myc-tag Magnetic Beads suspension to increase protein yield.
11. Incubate with gentle end-over-end mixing for 1 hour at 4°C.
12. Place the tube on Magnetic Rack (MBL; code no.3190) for 10 seconds. And remove the supernatant.
13. Add 1 mL of Wash Solution to the tube. Resuspend the magnetic bead pellet by vortexing for a few seconds.
14. Place the tube on the Magnetic Rack for 10 seconds. And remove the supernatant.
15. Repeat steps 13-14 three additional times.
16. Add 20 μ L Elution Peptide Solution to the tube. Vortex the tube for a few seconds. Incubate for 5 minutes at room temperature.
17. Place the tube on the Magnetic Rack for 10 seconds. Transfer the supernatant to another tube.
18. Repeat steps 16-17 again. The two eluates may be pooled in one sampling tube.

B. Purification from culture supernatant

(c-Myc-tagged protein is secreted into the culture supernatant)


1. Collect the culture supernatant from the cell culture dish into 15 mL centrifuge tube.
2. Centrifuge at 400 x g for 5 minutes to remove cell debris.
3. Transfer the supernatant to new 15 mL centrifuge tube.
4. Resuspend Anti-c-Myc-tag Magnetic Beads by vortexing for a few seconds immediately before dispensing.
5. Dispense 50 μ L of the Anti-c-Myc-tag Magnetic Beads suspension into the 15 mL centrifuge tube with culture supernatant. If necessary, please use 100 μ L of the Anti-c-Myc-tag Magnetic Beads suspension to increase protein yield
6. Incubate with gentle end-over-end mixing for 1 hour at 4°C.
7. Centrifuge at 400 x g for 5 minutes using swinging bucket rotor. Carefully remove the supernatant, but leave 100-400 μ L supernatant above the magnetic bead pellet.
8. Resuspend the magnetic bead pellet in the 100-400 μ L supernatant by vortexing for a few seconds.

- Transfer the suspension to 1.5 mL sampling tube.
9. Place the tube on Magnetic Rack (MBL; code no. 3190) for 10 seconds. Remove the supernatant.
 10. Add 1 mL of Wash Solution to the tube. Resuspend the magnetic bead pellet by vortexing for a few seconds.
 11. Place the tube on the Magnetic Rack for 10 seconds. Remove the supernatant.
 12. Repeat steps 10-11 three additional times.
 13. Add 20 μ L Elution Peptide Solution to the tube. Vortex the tube for a few seconds. Incubate for 5 minutes at room temperature.
 14. Place the tube on the Magnetic Rack for 10 seconds. Transfer the supernatant to another tube.
 15. Repeat steps 13-14 again. The two eluates may be pooled in one sampling tube.

(Note: Steps 14-18 of protocol A are identical to steps 11-15 of protocol B.)

Example of Purification Result 1

Purification of c-Myc-tagged GST protein

Lane 1: Input


Lane 2: Eluate (purified using c-Myc-tagged Protein Magnetic PURIFICATION KIT)

SDS-PAGE (Coomassie Brilliant Blue Staining)

The transformed competent *E.coli* JM109 cells were induced by the addition of IPTG in the culture medium (1 mL) to express the C-terminal c-Myc-tagged GST protein. After an induction period, the cells were pelleted by centrifugation and then resuspended in 0.5 mL of the Lysis buffer. c-Myc-tagged GST protein was purified according to the preceding protocol A.

Example of Purification Result 2

Purification of c-Myc-tagged Sox2 protein

Lane 1: Input

Lane 2: Eluate (purified using c-Myc-tagged Protein Magnetic PURIFICATION KIT)

SDS-PAGE (Coomassie Brilliant Blue Staining) & X-gal staining

Human embryonic kidney (HEK293T) cells were transfected with pcDNA-c-Myc-tagged Sox2 and cultured for 60 hours. Cells were then lysed in the Lysis buffer (0.5 mL/100-mm dish) and purified according to the preceding protocol A.

Additional Information

Several reagents were examined whether or not they were suitable for use with the c-Myc tagged Protein Magnetic PURIFICATION KIT. The results are listed below.

Chaotropic agents

	Urea	1 M	Yes
	Guanidine-HCl	1 M	Yes

Reducing agents

	DTT	10 mM	Yes
	2-Mercaptoethanol	10 mM	Yes

Surfactants

Nonionic	Tween-20	5%	Yes
	TritonX-100	5%	Yes
	NP-40	5%	Yes
	Digitonin	1%	Yes
	n-Octyl-β-D-gulcoside	1%	Yes
Zwitterionic	CHAPS	1%	Yes
	CHAPSO	1%	Yes
Anionic	SDS	0.05%	No
	Sodium Deoxycholate	0.5%	Yes

Others

	NaCl	1 M	Yes
	Glycerol	10%	Yes
	EDTA	10 mM	Yes

The “Yes” indicates the reagents can be used in the Lysis buffer for this kit up to the indicated concentration. The “No” indicates the reagents cannot be used in the Lysis buffer for this kit at the indicated concentration.

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