

ELISA Kit for Measuring Human soluble Nectin-4

CircuLex Human soluble Nectin-4 ELISA Kit

Cat# CY-8217

Intended Use.....	1
Storage.....	1
Introduction.....	2
Principle of the Assay.....	3
Materials Provided.....	4
Materials Required but not Provided.....	4
Precautions and Recommendations.....	5
Sample Collection and Storage.....	5
Detailed Protocol.....	6-8
Calculations.....	9
Measurement Range.....	9
Troubleshooting.....	9
Reagent Stability.....	9
Assay Characteristics.....	10-12
Example of Test Results.....	13
References.....	14

Intended Use

The MBL Research Product **CircuLex Human soluble Nectin-4 ELISA Kit** is used for the quantitative measurement of the soluble form of human nectin-4 in cell culture supernatant.

Individual users should determine appropriate conditions when using other types of samples.

This assay kit is for research use only and not for use in diagnostic or therapeutic procedures.

Storage

- Upon receipt store all components at 4°C.
- Do not expose reagents to excessive light.

Introduction


Nectins are a family of type I transmembrane proteins that are involved in a Ca^{2+} -independent immunoglobulin-like cell-cell adhesion molecules consisting of four members, Nectin-1-4 (1). Nectin-4, also known as PVRL4, LNIR or PRR4, is expressed in the embryo and the placenta, whereas the other nectins are expressed widely in human adult tissues. The extracellular domain of nectin-4 is cleaved at the cell surface via the proteolytic activity of ADAM17/TACE, and the soluble form is released into the blood. Increased levels of the circulating form of nectin-4 have been detected in the sera of patients with metastatic breast cancer (2). The soluble form of nectin-4 has also been detected in significant amounts in the sera of patients with non-small cell lung cancer (NSCLC), and the up-regulation of soluble nectin-4 in the serum has been suggested to be associated with a poor prognosis for NSCLC (3).

Nectin-4 is overexpressed in a variety of cancers (4-6) other than the above and is currently under clinical investigation as a therapeutic target. Oncolytic virotherapy is a novel treatment strategy for cancer. Vaccine and laboratory-adapted strains of measles virus (MV) uses nectin-4 and CD46 molecule as major entry receptors into cells (7). A recombinant measles virus which selectively uses nectin-4 receptor is used in oncolytic agents for cancer (4,5,8). Also anti-nectin-4 antibodies such as enfortumab (AGS-22M6E) are under clinical investigation for cancer therapies (9,10). Furthermore mutations in nectin-4 are associated with ectodermal dysplasia-syndactyly syndrome 1 (EDSS1) (11,12).

Principle of the Assay

The MBL Research Product **CircuLex Human soluble Nectin-4 ELISA Kit** employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human nectin-4 is pre-coated onto a microplate. Standards and samples are pipetted into the wells and the immobilized antibody binds any human soluble nectin-4 present. After washing away any unbound substances, a biotinylated monoclonal antibody specific for human nectin-4 is added to the wells. Following a wash to remove any unbound biotinylated antibody, Streptavidin-HRP is added to the wells. After washing away any unbound Streptavidin-HRP, HRP remaining on the well is allowed to react with the substrate H₂O₂-tetramethylbenzidine. The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured at 450 nm. The absorbance is proportional to the concentration of human soluble nectin-4. A standard curve is constructed by plotting absorbance values versus human soluble nectin-4 concentrations of calibrators, and concentrations of unknown samples are determined using this standard curve.

Summary of Procedure



Materials Provided

All samples and standards should be assayed in duplicate. The following components are supplied and are sufficient for the one 96-well microplate kit.

Microplate: One microplate supplied ready to use, with 96 wells (12 strips of 8-wells) in a foil, zip-lock bag with a desiccant pack. Wells are pre-coated with anti-human nectin-4 monoclonal antibody as a capture antibody.

10X Wash Buffer: One bottle containing 100 mL of 10X buffer containing Tween®-20.

Sample Dilution Buffer: One bottle containing 50 mL of 1X buffer; use for standard and sample dilution. Ready to use.

Human soluble Nectin-4 Standard: One vial containing X* ng of lyophilized recombinant human soluble nectin-4.

***The amount is changed depending on lot. See the real "User's Manual" included in the kit box.**

Biotinylated Detection Antibody: One vial containing 12 mL of biotinylated anti-human nectin-4 monoclonal antibody. Ready to use.

100X Streptavidin-HRP: One vial containing 120 µL of 100X HRP (horseradish peroxidase) conjugated streptavidin.

Streptavidin-HRP Dilution Buffer: One bottle containing 12 mL of 1X buffer; use for dilution of 100X Streptavidin-HRP.

Substrate Reagent: One bottle containing 20 mL of the chromogenic substrate, tetra-methylbenzidine (TMB). Ready to use.

Stop Solution: One bottle containing 20 mL of 1 N H₂SO₄. Ready to use.

Materials Required but not Provided

- **Pipettors:** 2-20 µL, 20-200 µL and 200-1,000 µL precision pipettors with disposable tips.
- **Precision repeating pipettor**
- **Orbital microplate shaker**
- **Microcentrifuge and tubes** for sample preparation.
- **Vortex mixer**
- **(Optional) Microplate washer:** Manual washing is possible but not preferable.
- **Plate reader** capable of measuring absorbance in 96-well plates at dual wavelengths of 450 nm/540 nm. Dual wavelengths of 450/550 or 450/595 nm can also be used. The plate can also be read at a single wavelength of 450 nm, which will give a somewhat higher reading.
- **(Optional) Software package facilitating data generation and analysis**
- **500 or 1,000 mL graduated cylinder.**
- **Reagent reservoirs**
- **Deionized water of the highest quality**
- **Disposable paper towels**

Precautions and Recommendations

- **Although we suggest to conduct experiments as outlined below, the optimal experimental conditions will vary depending on the parameters being investigated, and must be determined by the individual user.**
- Allow all the components to come to room temperature before use.
- All microplate strips that are not immediately required should be returned to the zip-lock pouch, which must be carefully resealed to avoid moisture absorption.
- Do not use kit components beyond the indicated kit expiration date.
- Use only the microtiter wells provided with the kit.
- Rinse all detergent residue from glassware.
- Use deionized water of the highest quality.
- Do not mix reagents from different kits.
- The buffers and reagents in this kit may contain preservatives or other chemicals. Care should be taken to avoid direct contact with these reagents.
- Do not mouth pipette or ingest any of the reagents.
- Do not smoke, eat, or drink when performing the assay or in areas where samples or reagents are handled.
- Dispose of tetra-methylbenzidine (TMB) containing solutions in compliance with local regulations.
- Avoid contact with Substrate Solution which contains hydrogen peroxide.
- **CAUTION: Biological samples may be contaminated with infectious agents. Do not ingest, expose to open wounds or breathe aerosols. Wear protective gloves and dispose of biological samples properly.**
- **CAUTION: Stop Solution is a strong acid. Wear disposable gloves and eye protection when handling the solution.**

Sample Collection and Storage

Cell culture supernatant: Remove any particulates by centrifugation and assay immediately or aliquot and store samples at below -70°C. Avoid repeated freeze-thaw cycles.

Detailed Protocol

The MBL Research Product **CircuLex Human soluble Nectin-4 ELISA Kit** is provided with removable strips of wells so the assay can be carried out on separate occasions using only the number of strips required for the particular determination. Since experimental conditions may vary, an aliquot of the standard within the kit, should be included in each assay as a calibrator. Disposable pipette tips and reagent troughs should be used for all liquid transfers to avoid cross-contamination of reagents or samples.

Preparation of Working Solutions

All reagents need to be brought to room temperature prior to the assay. Assay reagents are supplied ready-to-use, with the exception of **10X Wash Buffer**, **100X Streptavidin-HRP**, and **Human soluble Nectin-4 Standard**.

1. Prepare a working solution of **Wash Buffer** by adding **100 mL** of **10X Wash Buffer** to **900 mL** of deionized (distilled) water (**ddH₂O**). Mix well. Store at 4°C for two weeks or -20°C for long-term storage.
2. Prepare **1X Streptavidin-HRP** by 100-fold diluting **100X Streptavidin-HRP** with **Streptavidin-HRP Dilution Buffer** at the time of use and discard any unused portion after use.
3. Reconstitute **Human soluble Nectin-4 Standard** with **X* μL** of **ddH₂O** by gently mixing. After reconstitution, immediately dispense it in small aliquots (e.g. 100 μL) to plastic micro-centrifuge tubes and store below -70°C to avoid non-specific adsorption to glass surface and multiple freeze-thaw cycles. The concentration of human soluble nectin-4 in vial should be **19.2 ng/mL**, which is referred to as the **Master Standard** of human soluble nectin-4.

***The amount is changed depending on lot. See the real "User's Manual" included in the kit box.**

Prepare Standard Solutions as follows:

Use **Master Standard** to produce a dilution series (below). Mix each tube thoroughly before the next transfer. **Std.1 (1,920 pg/mL)** serves as the highest standard. **Sample Dilution Buffer** serves as the zero standard (Blank).

	Volume of Standard	Sample Dilution Buffer	Concentration
Std.1	60 μL of Master Standard (19.2 ng/mL)	540 μL	1,920 pg/mL
Std.2	300 μL of Std. 1 (1,920 pg/mL)	300 μL	960 pg/mL
Std.3	300 μL of Std. 2 (960 pg/mL)	300 μL	480 pg/mL
Std.4	300 μL of Std. 3 (480 pg/mL)	300 μL	240 pg/mL
Std.5	300 μL of Std. 4 (240 pg/mL)	300 μL	120 pg/mL
Std.6	300 μL of Std. 5 (120 pg/mL)	300 μL	60 pg/mL
Std.7	300 μL of Std. 6 (60 pg/mL)	300 μL	30 pg/mL
Blank	-	300 μL	0 pg/mL

Note: Do not use a repeating pipette. Change tips for every dilution. Wet tip with Dilution Buffer before dispensing.

Sample Preparation

Dilute samples with **Sample Dilution Buffer**.

- Cell culture supernatants require appropriate dilutions.

Assay Procedure

1. Remove the appropriate number of microtiter wells from the foil pouch and place them into the well holder. Return any unused wells to the foil pouch, refold, seal with tape and store at 4°C.
2. Dilute samples with **Sample Dilution Buffer**. (See “Sample Preparation” above.)
3. Pipette **100 µL** of **Standard Solutions (Std1-Std7, Blank)** and **diluted samples** in duplicates, into the appropriate wells.
4. Incubate the plate **at room temperature (ca.25°C) for 60 minutes**, shaking at ca. 300 rpm on an orbital microplate shaker.
5. Wash 4-times by filling each well with Wash Buffer (350 µL) using a squirt bottle, multi-channel pipette, manifold dispenser or microplate washer.
6. Add **100 µL** of **Biotinylated Detection Antibody** into each well.
7. Incubate the plate **at room temperature (ca.25°C) for 60 minutes**, shaking at ca. 300 rpm on an orbital microplate shaker.
8. Wash 4-times by filling each well with Wash Buffer (350 µL) using a squirt bottle, multi-channel pipette, manifold dispenser or microplate washer.
9. Add **100 µL** of **1X Streptavidin-HRP** into each well.
10. Incubate the plate **at room temperature (ca.25°C) for 20 minutes**, shaking at ca. 300 rpm on an orbital microplate shaker.
11. Wash 4-times by filling each well with Wash Buffer (350 µL) using a squirt bottle, multi-channel pipette, manifold dispenser or microplate washer.
12. Add **100 µL** of **Substrate Reagent**. Avoid exposing the microtiter plate to direct sunlight. Covering the plate with e.g. aluminum foil is recommended. Return Substrate Reagent to 4°C immediately after the necessary volume is removed
13. Incubate the plate **at room temperature (ca.25°C) for 10-20 minutes**, shaking at ca. 300 rpm on an orbital microplate shaker. The incubation time may be extended up to 30 minutes if the reaction temperature is below than 20°C.
14. Add **100 µL** of **Stop Solution** to each well in the same order as the previously added Substrate Reagent.
15. Measure absorbance in each well using a spectrophotometric microplate reader at dual wavelengths of 450/540 nm. Dual wavelengths of 450/550 or 450/595 nm can also be used. Read the microplate at

450 nm if only a single wavelength can be used. Wells must be read within 30 minutes of adding the Stop Solution.

Note-1: Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.

Note-2: Reliable standard curves are obtained when either O.D. values do not exceed 0.25 units for the blank (zero concentration), or 3.0 units for the highest standard concentration.

Note-3: If the microplate reader is not capable of reading absorbance greater than the absorbance of the highest standard, perform a second reading at 405 nm. A new standard curve, constructed using the values measured at 405 nm, is used to determine the concentration of off-scale samples. The readings at 405 nm should not replace the on-scale readings at 450 nm.

Calculations

Average the duplicate readings for each standard, control and sample, and subtract the optical density of the average zero standard. Plot the optical density versus the concentration of standards and draw the best curve. Most microtiter plate readers perform automatic calculations of analyte concentration. The standard curve fits best to a sigmoidal four-parameter logistic equation. The results of unknown samples can be calculated with any computer program having a four-parameter logistic function.

A standard curve is also to be constructed by plotting the absorbance (Y) versus log of the known concentration (X) of standards, using a cubic function. Alternatively, the logit log function can be used to linearize the standard curve (i.e. logit of optical density (Y) is plotted versus log of the known concentration (X) of standards). To determine the concentration of each sample, first find the optical density on the y-axis and extend a horizontal line to the standard curve. At the point of intersection, extend a vertical line to the x-axis and read the corresponding concentration.

If the samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

Measurement Range

The measurement range is 30 pg/mL to 1,920 pg/mL. Any sample reading higher than the highest standard should be diluted with Sample Dilution Buffer in higher dilution and re-assayed. Dilution factors need to be taken into consideration in calculating the human soluble nectin-4 concentration.

Troubleshooting

1. All samples and standards should be assayed in duplicate, using the protocol described in **Detailed Protocol**. Incubation times or temperatures significantly different from those specified may give erroneous results.
2. Poor duplicates, accompanied by elevated values for wells containing no sample, indicate insufficient washing. If all instructions in **Detailed Protocol** were followed accurately, such results indicate a need for washer maintenance.
3. Overall low signal may indicate that desiccation of the plate has occurred between the final wash and addition of Substrate Reagent. Do not allow the plate to dry out. Add Substrate Reagent immediately after wash.

Reagent Stability

All of the reagents included in the MBL Research Product **CircuLex Human soluble Nectin-4 ELISA Kit** have been tested for stability. Reagents should not be used beyond the stated expiration date.


Assay Characteristics

1. Sensitivity

The limit of detection (defined as such a concentration of human soluble nectin-4 giving absorbance higher than mean absorbance of blank* plus three standard deviations of the absorbance of blank: A blank + 3SD blank) is better than 5.8 pg/mL of sample.

* Sample Dilution Buffer was pipetted into blank wells.

Typical standard curve



2. Precision

Intra-assay Precision (Precision within an assay)

Three samples* of known concentration were tested sixteen times on one plate to assess intra-assay precision.

- Intra-assay (Within-Run, n=16) CV=3.7, 3.0, 3.8 %

* Sample: Cell culture supernatant

Human nectin-4 conc. (ng/mL)			
	Sample 1	Sample 2	Sample 3
1	0.35	0.86	1.37
2	0.35	0.81	1.27
3	0.34	0.81	1.26
4	0.35	0.89	1.33
5	0.35	0.83	1.39
6	0.34	0.84	1.30
7	0.35	0.78	1.30
8	0.36	0.83	1.41
9	0.37	0.83	1.41
10	0.37	0.83	1.37
11	0.37	0.82	1.39
12	0.38	0.85	1.40
13	0.36	0.84	1.40
14	0.37	0.85	1.39
15	0.37	0.87	1.41
16	0.38	0.83	1.36
MAX.	0.38	0.89	1.41
MIN.	0.34	0.78	1.26
MEAN	0.36	0.84	1.36
S.D.	0.01	0.03	0.05
C.V.	3.7%	3.0%	3.8%

Inter-assay Precision (Precision between assays)

Three samples* of known concentration were tested in five separate assays to assess inter-assay precision.

- Inter-assay (Run-to-Run, n=5) CV=6.5, 4.2, 6.5 %


* Sample: Cell culture supernatant

Human nectin-4 conc. (ng/mL)			
	Serum 1	Serum 2	Serum 3
1	2.34	6.03	13.95
2	2.30	5.67	14.35
3	2.09	5.65	13.22
4	2.04	5.56	13.83
5	2.33	6.09	14.23
MAX.	2.34	6.09	14.35
MIN.	2.04	5.56	13.22
MEAN	2.22	5.80	13.92
S.D.	0.144	0.242	0.441
C.V.	6.5%	4.2%	3.2%

3. Linearity


To assess the linearity of the assay, three samples* were serially diluted with Sample Dilution Buffer to produce samples with values within the dynamic range of the assay.

* Sample: Cell culture supernatant



Example of Test Results

Fig. Human soluble nectin-4 concentration in cell culture supernatants



References

1. Nakanishi H and Takai Y. Roles of nectins in cell adhesion, migration and polarization. *Biol Chem.* 2004 Oct;385(10):885-92.
2. Fabre-Lafay S *et al.* Nectin-4, a new serological breast cancer marker, is a substrate for tumor necrosis factor-alpha-converting enzyme (TACE)/ADAM-17. *J Biol Chem.* 2005 May 20;280(20):19543-50.
3. Takano A *et al.* Identification of nectin-4 oncoprotein as a diagnostic and therapeutic target for lung cancer. *Cancer Res.* 2009 Aug 15;69(16):6694-703.
4. Amagai Y *et al.* Oncolytic Activity of a Recombinant Measles Virus, Blind to Signaling Lymphocyte Activation Molecule, Against Colorectal Cancer Cells. *Sci Rep.* 2016 Apr 19;6:24572.
5. Fujiyuki T *et al.* A measles virus selectively blind to signaling lymphocytic activation molecule shows anti-tumor activity against lung cancer cells. *Oncotarget.* 2015 Sep 22;6(28):24895-903.
6. Nishiwada S *et al.* Nectin-4 expression contributes to tumor proliferation, angiogenesis and patient prognosis in human pancreatic cancer. *J Exp Clin Cancer Res.* 2015 Mar 28;34:30.
7. Noyce RS and Richardson CD. Nectin 4 is the epithelial cell receptor for measles virus. *Trends Microbiol.* 2012 Sep;20(9):429-39.
8. Shoji K *et al.* Development of new therapy for canine mammary cancer with recombinant measles virus. *Mol Ther Oncolytics.* 2016 Jan 13;3:15022.
9. Challita-Eid PM *et al.* Enfortumab Vedotin Antibody-Drug Conjugate Targeting Nectin-4 Is a Highly Potent Therapeutic Agent in Multiple Preclinical Cancer Models. *Cancer Res.* 2016 May 15;76(10):3003-13.
10. Campbell DO *et al.* Preclinical Evaluation of an Anti-Nectin-4 ImmunoPET Reagent in Tumor-Bearing Mice and Biodistribution Studies in Cynomolgus Monkeys. *Mol Imaging Biol.* 2016 Oct;18(5):768-75.
11. Brancati F *et al.* Mutations in PVRL4, encoding cell adhesion molecule nectin-4, cause ectodermal dysplasia-syndactyly syndrome. *Am J Hum Genet.* 2010 Aug 13;87(2):265-73.
12. Raza SI *et al.* A novel homozygous nonsense mutation in the PVRL4 gene and expansion of clinical spectrum of EDSS1. *Ann Hum Genet.* 2015 Mar;79(2):92-8.

For more information, please visit our web site.

<https://ruo.mbl.co.jp/>

MANUFACTURED BY

MBL A JSR Life Sciences Company
MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.
URL: <https://ruo.mbl.co.jp>
E-mail: support@mbl.co.jp

CycLex/CircuLex products are supplied for research use only. CycLex/CircuLex products and components thereof may not be resold, modified for resale, or used to manufacture commercial products without prior written approval from MBL. To inquire about licensing for such commercial use, please contact us via email.