



cGK Positive Control (Catalytic domain)

Product Data Sheet

For Research Use Only, Not for use in diagnostic procedures



cGK Positive Control (Catalytic domain) **(Human, a.a.324-671, recombinant protein expressed in Sf9 cells)** **Cat# CY-E1161-1**

Lot No.
For 200 Assays
400 units (4 units/ μ L)

Product Description: The catalytic domain of human cGK (cGMP-Dependent Protein Kinase), containing an N-terminal GST tag and a C-terminal His tag, is expressed in recombinant baculovirus infected Sf9 cells and purified by sequentially using GSH agarose and Ni-NTA agarose chromatography. The cGK Positive Control (Catalytic domain) is designed to use for CycLex Cyclic GMP dependent protein kinase (cGK) Assay Kit (Cat# CY-1161). The cGK Positive Control (Catalytic domain) should be added to the well at 2 units/well. This enzyme shows protein kinase activity in the absence of cGMP under the protocol described below. Unused cGK Positive Control (Catalytic domain) should be stored at -70°C .

Product Size: 400 units/100 μ L

Formulation: Supplied frozen in a buffer containing 20 mM Hepes-KOH (pH 7.5), 1 % BSA, 1 mM EDTA, 1 mM DTT, 50 mM NaCl, 0.03 % Brij35 and 50 % glycerol.

Source: Human cGMP-Dependent Protein Kinase, 324-671, containing an N-terminal GST tag and a C-terminal His tag, expressed in Sf9 cells.

Molecular Weight: 64 kDa band by SDS-PAGE analysis.

Purity: > 90 % pure as determined by SDS-PAGE analysis.

Substrates: cGMP-Dependent Protein Kinase phosphorylates a number of substrates, including histone proteins H2b and H1, brain G protein and high mobility group 14 protein.

Inhibitors: Polycations, such as poly-L-arginine, inhibit cGMP-Dependent Protein Kinase.

Unit Definition: One unit is defined as the amount of kinase required to incorporate 1pmol of phosphate into the GST-G substrate fusion protein, per minute at 30°C .

Assay Conditions: Assay activity of cGMP-Dependent Protein Kinase in a 50 μ L reaction containing 20 mM Hepes KOH (pH 7.5), 5 mM MgCl_2 , 1 mM DTT, 100 μ M [γ ^{32}P] ATP (1 μ Ci) and 4 μ g of GST-G substrate fusion protein in the absence of cGMP. Start the reaction by adding 10 μ L of the enzyme, diluted 50-fold in a buffer containing 20 mM Hepes KOH (pH 7.5), 1 mM DTT, 0.03 % Brij35. Incubate for 30 minutes at 30°C . Terminate the reaction by adding 600 μ L of cold 10 % TCA solution containing 0.2 % sodium pyrophosphate and stand on ice for 15 min. Filtrate acid insoluble material through GFC filters (Whatman Inc.), wash 4 times with 1 % TCA and rinse filters with ethanol. Dry filters and count in a liquid scintillation counter.



cGK Positive Control (Catalytic domain)

Product Data Sheet

For Research Use Only, Not for use in diagnostic procedures



Storage Conditions: Stable for 12 months at -70°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot enzyme to avoid repeated freezing and thawing.

References:

1. Edelman, A.M., Blumenthal, D.K. and Krebs, E.G. Protein serine/threonine kinases. *Ann. Rev. Biochem.* 56, 567, 1987
2. Beebe, S.J. and Corbin, J.D. In: *The Enzymes*, Vol. 17, 3rd ed., Boyer, P.D. and Krebs, E.D., eds., 44, 1986
3. Corbin, J.D. and Doskeland, S.O. *J. Biol. Chem.* 258, 11391, 1983
4. Endo, S et al. *Proc. Natl. Acad. Sci. USA.* 96, 2467-2472, 1999

For more information, please visit our web site.

<https://ruo.mbl.co.jp/>

MANUFACTURED BY

MBL A JSR Life Sciences Company

MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.

URL: <https://ruo.mbl.co.jp>

E-mail: support@mbl.co.jp

CycLex/CircuLex products are supplied for research use only. CycLex/CircuLex products and components thereof may not be resold, modified for resale, or used to manufacture commercial products without prior written approval from MBL. To inquire about licensing for such commercial use, please contact us via email.