

POLYCLONAL ANTIBODY

# Anti- $\beta$ -galactosidase pAb

Code No.	Quantity	Form
PM049	100 $\mu$ L	Rabbit IgG

**BACKGROUND:**  $\beta$ -galactosidase is a homo-tetrameric enzyme, with each subunit having a molecular weight of 116 kDa. Eukaryotic genes are often expressed as fusion protein by the  $\beta$ -galactosidase (*lacZ*) gene, resulting in the expression of a fusion hybrid with  $\beta$ -galactosidase. Anti- $\beta$ -galactosidase antibody provides a simple method to isolate fusion proteins directly from crude bacterial lysates, using immunoaffinity chromatography or immunoprecipitation. Anti- $\beta$ -galactosidase can also be used for the immunocytochemical detection of  $\beta$ -galactosidase in cells and tissues that express transfected bacterial *lacZ* gene or  $\beta$ -galactosidase fusion protein.

**SOURCE:** This antibody was purified from rabbit serum using protein A agarose. The rabbit was immunized with full length *E. coli*  $\beta$ -galactosidase.

**FORMULATION:** 100  $\mu$ L volume of PBS containing 50% glycerol, pH 7.2. No preservative is contained.

**STORAGE:** This antibody solution is stable for one year from the date of purchase when stored at  $-20^{\circ}\text{C}$ .

**REACTIVITY:** This antibody reacts with  $\beta$ -galactosidase on Western blotting, Immunoprecipitation, Immunohistochemistry and Immunocytochemistry.

## APPLICATIONS:

Western blotting: 1:1,000

Immunoprecipitation: 1  $\mu$ L/sample

Immunohistochemistry: 1:200

Immunocytochemistry: 1:100

Flow cytometry: Not tested

Detailed procedure is provided in the following **PROTOCOLS**.

## INTENDED USE:

For Research Use Only. Not for use in diagnostic procedures.

## REFERENCES:

- 1) Otsu, K., *et al.*, *J. Bone Miner. Res.* **31**, 1943-1954 (2016) [IHC]
- 2) Vaish, V., *et al.*, *Genes Chromosomes Cancer.* **55**, 577-590 (2016)




**Western blot analysis of  $\beta$ -galactosidase expression in transfectant using PM049.**

The descriptions of the following protocols are examples. Each user should determine the appropriate condition.

## PROTOCOLS:

### SDS-PAGE & Western Blotting

- 1) Wash the  $1 \times 10^6$  transfectant cells 3 times with PBS and suspend with 1 mL of Laemmli's sample buffer.
- 2) Boil the samples for 2 minutes and centrifuge. Load 10  $\mu$ L of the sample per lane in a 1 mm thick SDS-polyacrylamide gel for electrophoresis.
- 3) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at 1 mA/cm<sup>2</sup> for 1 hour in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacturer's manual for precise transfer procedure.
- 4) To reduce nonspecific binding, soak the membrane in 10% skimmed milk (in PBS, pH 7.2) for 1 hour at room temperature, or overnight at  $4^{\circ}\text{C}$ .
- 5) Incubate the membrane with primary antibody diluted with PBS, pH 7.2 containing 1% skimmed milk as suggested in the **APPLICATIONS** for 1 hour at room temperature. (The concentration of antibody will depend on condition.)
- 6) Wash the membrane with PBS-T [0.05% Tween-20 in PBS] (5 minutes x 3 times).
- 7) Incubate the membrane with 1:10,000 of Anti-IgG (Rabbit) pAb-HRP (MBL; code no. 458) diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hour at room temperature.
- 8) Wash the membrane with PBS-T (5 minutes x 3 times).
- 9) Wipe excess buffer on the membrane, then incubate it with appropriate chemiluminescence reagent for 1 minute. Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
- 10) Expose to an X-ray film in a dark room for 3 minutes. Develop the film as usual. The condition for exposure and development may vary.



**Immunoprecipitation of  $\beta$ -galactosidase from transfectant with normal rabbit IgG (1) or PM049 (2). After immunoprecipitated with the antibody, immunocomplex was resolved on SDS-PAGE and immunoblotted with anti- $\beta$ -galactosidase monoclonal antibody (MBL; code no. M094-3).**


#### **Immunoprecipitation**

- 1) Wash the  $5 \times 10^6$  transfectant cells 3 times with PBS and suspend with 1 mL of cold Lysis buffer [50 mM Tris-HCl (pH 7.5), 150 mM NaCl, 0.05% NP-40] containing appropriate protease inhibitors. Incubate it at 4°C with rotating for 15 minutes, then sonicate briefly (up to 10 seconds).
- 2) Centrifuge the tube at 12,000 x g for 10 minutes at 4°C and transfer the supernatant to another tube.
- 3) Add primary antibody as suggested in the **APPLICATIONS** into 200  $\mu$ L of the supernatant. Mix well and incubate with gentle agitation for 60-120 minutes at 4°C.
- 4) Add 20  $\mu$ L of 50% protein A agarose beads resuspended in the cold Lysis buffer. Mix well and incubate with gentle agitation for 60 minutes at 4°C.
- 5) Wash the beads 3-5 times with the cold Lysis buffer (centrifuge the tube at 2,500 x g for 10 seconds).
- 6) Resuspend the beads in 20  $\mu$ L of Laemmli's sample buffer, boil for 3-5 minutes, and centrifuge for 5 minutes. Use 10  $\mu$ L/lane for the SDS-PAGE analysis.  
(See **SDS-PAGE & Western blotting**.)

#### **Immunohistochemical staining for paraffin-embedded sections**

- 1) Deparaffinize the sections with Xylene 3 times for 3-5 minutes each.
- 2) Wash the slides with Ethanol 3 times for 3-5 minutes each.
- 3) Wash the slides with PBS 3 times for 3-5 minutes each.
- 4) Remove the slides from the PBS and cover each section with 3% H<sub>2</sub>O<sub>2</sub> for 10 minutes at room temperature to block endogenous peroxidase activity. Wash 3 times in PBS for 5 minutes each.
- 5) Remove the slides from PBS, wipe gently around each section and cover tissues with blocking buffer (20 mM HEPES, 1% BSA, 135 mM NaCl) for 5 minutes to block non-specific staining. Do not wash.

- 6) Tip off the blocking buffer, wipe gently around each section and cover tissues with primary antibody diluted with blocking buffer as suggested in the **APPLICATIONS**.
- 7) Incubate the sections for 1 hour at room temperature.
- 8) Wash the slides 3 times in PBS for 5 minutes each.
- 9) Wipe gently around each section and cover tissues with ENVISION+Dual Link (DAKO; code no. K4063). Incubate for 1 hour at room temperature. Wash as in step 9).
- 10) Visualize by reacting for 10 minutes with DAB substrate solution (DAKO; code no. K3465). \*DAB is a suspect carcinogen and must be handled with care. Always wear gloves.
- 11) Wash the slides in water for 5 minutes.
- 12) Counter stain in hematoxylin for 1 minute, wash the slides 3 times in water for 5 minutes each, and then immerse the slides in PBS for 5 minutes. Dehydrate by immersing in Ethanol 3 times for 3 minutes each, followed by immersing in Xylene 3 times for 3 minutes each.
- 13) Now ready for mounting.




#### ***Immunohistochemical detection of $\beta$ -galactosidase in frozen section of Lewis lung carcinoma xenograft with PM049.***

This data was provided by Dr. Minami (Laboratory for Systems Biology and Medicine at RCAST, University of Tokyo)


#### **Immunohistochemical staining for frozen sections: For acetone fixed section**

- 1) Wash the slide in PBS (5 minutes x 3 times).
- 2) Immerse the slide in Image iTTM FX™ for 30 minutes at room temperature
- 3) Wash the slide with PBS 3 times for 5 minutes each.
- 4) Remove the slides from PBS, wipe gently around each section and cover tissues with Protein Block (Dako; code no. X0909) for 20 minutes to block non-specific staining. Do not wash.
- 5) Tipp off the blocking buffer, wipe gently around each section and cover tissues with primary antibody diluted with Protein Block as suggested in the **APPLICATIONS**.
- 6) Incubate the sections overnight at 4°C.
- 7) Wash the slide with PBS 3 times for 5 minutes each.

- 8) Wipe gently around each section and cover tissues with 1:50 Alexa Fluor® 488 conjugated anti-rabbit IgG (Invitrogen; code no. A11008). Incubate for 1 hour at room temperature.
- 9) Wash the slide with PBS 3 times for 5 minutes each.
- 10) Promptly add mounting medium onto the slide, then put a cover slip on it.
- 10) The glass slide was washed 3 times with PBS.
- 11) Wipe excess liquid from slide but take care not to touch the cells. Never leave the cells to dry.
- 12) Promptly add mounting medium onto the slide, then put a cover slip on it.



***Immunocytochemical detection of  $\beta$ -galactosidase in transfectant with PM049.***



***Immunocytochemical detection of  $\beta$ -galactosidase in paraffin embedded section of transfectant with PM049.***

#### **Immunocytochemistry**

- 1) Culture the cells in the appropriate condition on a glass slide. (for example, spread  $1 \times 10^4$  cells for one slide, then incubate in a CO<sub>2</sub> incubator for one night.)
- 2) Wash the cells 3 times with PBS.
- 3) Fix the cells by immersing the slide in PBS containing 4% paraformaldehyde for 10 minutes at room temperature.
- 4) The glass slide was washed with PBS 3 times.
- 5) Immerse the slide in PBS containing 0.1% TritonX-100 for 10 minutes at room temperature.
- 6) The glass slide was washed 3 times with PBS.
- 7) Add the primary antibody diluted with PBS as suggested in the **APPLICATIONS** onto the cells and incubate for 30 minutes at room temperature (Optimization of antibody concentration or incubation condition are recommended if necessary.)
- 8) The glass slide was washed 3 times with PBS.
- 9) Add 100  $\mu$ L of 1:500 Alexa Fluor® 488 conjugated anti-rabbit IgG (Invitrogen; code no. A11008) diluted with PBS onto the cells. Incubate for 30 minutes at room temperature. Keep out light by aluminum foil.