

# Normal Guinea Pig IgG

<b>CODE No.</b>	PM067
<b>CLONALITY</b>	Polyclonal
<b>QUANTITY</b>	100 µL, 1 mg/mL
<b>SOURCE</b>	Purified IgG from normal guinea pig serum using protein A agarose.
<b>REACTIVITY</b>	No specific reaction was detected on immunoprecipitation and flow cytometry.
<b>FORMULATION</b>	1 mg/mL in PBS containing 50% glycerol. No preservative is contained.
<b>STORAGE</b>	This antibody solution is stable for one year from the date of purchase when stored at -20°C.

## APPLICATIONS-CONFIRMED

### Immunoprecipitation

### Flow cytometry

This antibody can be used as a negative isotypic control.  
The concentration will depend on the conditions.

## APPLICATION-REPORTED

Immunohistochemistry Reference 1) and 2)

<b>REFERENCES</b>	1) Murai, N., <i>et al.</i> , <i>PLoS One</i> <b>12</b> , e0186637 (2017) [IHC] 2) Yamane, T., <i>et al.</i> , <i>PLoS One</i> <b>12</b> , e0176809 (2017) [WB, IHC]
-------------------	---

For more information, please visit our web site <https://ruo.mbl.co.jp/>

## RELATED PRODUCTS

Please visit our web site <https://ruo.mbl.co.jp/>

The descriptions of the following protocols are examples. Each user should determine the appropriate condition.

### **Immunoprecipitation**

- 1) Wash  $1 \times 10^7$  cells 2 times with PBS and resuspend them with 1 mL of ice-cold Lysis buffer [50 mM Tris-HCl (pH 7.4), 150 mM NaCl, 0.05% NP-40] containing appropriate protease inhibitors. Incubate it on ice for 15 min., thereafter, sonicate briefly (up to 15 sec.).
- 2) Centrifuge the tube at  $12,000 \times g$  for 5 min. at  $4^\circ\text{C}$  and transfer the supernatant to another tube.
- 3) Add the isotype control antibody at the equal amount of the antibody for immunoprecipitation to the supernatant. Vortex briefly and incubate with gentle agitation for 60-120 min. at  $4^\circ\text{C}$ .
- 4) Mix 20  $\mu\text{L}$  of 50% protein A agarose beads slurry resuspended in 400  $\mu\text{L}$  of IP buffer [10 mM Tris-HCl (pH 8.0), 500 mM NaCl, 0.1% NP-40] with primary antibody. Incubate with gentle agitation for 1 hr. at room temperature.
- 5) Wash the beads 3 times with 1 mL of IP buffer.
- 6) Add 300  $\mu\text{L}$  of cell lysate (prepared sample of step 2), then incubate with gentle agitation for 1 hr. at room temperature.
- 7) Centrifuge the tube at  $2,500 \times g$  for 10 sec. Carefully remove and discard the supernatant.
- 8) Resuspend the beads with 1 mL of Lysis buffer.
- 9) Centrifuge the tube at  $2,500 \times g$  for 10 sec. Carefully remove and discard the supernatant.
- 10) Repeat Steps 8)-9) 5 times.
- 11) Resuspend the beads in 20  $\mu\text{L}$  of Laemmli's sample buffer, boil for 3 min. and centrifuge for 5 min.
- 12) Load 10  $\mu\text{L}$  of the sample per lane in a 1-mm-thick SDS-polyacrylamide gel for electrophoresis.
- 13) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at  $1 \text{ mA}/\text{cm}^2$  for 1 hr. in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacturer's manual for precise transfer procedure.
- 14) To reduce nonspecific binding, soak the membrane in 10% skimmed milk (in PBS, pH 7.2) for 1 hr. at room temperature.
- 15) Incubate the membrane with primary antibody diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hr. at room temperature. (The concentration of antibody will depend on the conditions.)
- 16) Wash the membrane with PBS (5 min.  $\times$  3 times).
- 17) Incubate the membrane with Anti-IgG (Rabbit) pAb-HRP (MBL; code no. 458) diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hr. at room temperature.
- 18) Wash the membrane with PBS (5 min.  $\times$  3 times).
- 19) Wipe excess buffer on the membrane, and then incubate it with appropriate chemiluminescence reagent for 1 min.
- 20) Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
- 21) Expose to an X-ray film in a dark room for 30 sec. Develop the film as usual. The condition for exposure and development may vary.



### ***Immunoprecipitation of mouse p62 from NIH/3T3 cells***


Lane 1: IP with Normal Guinea Pig IgG (PM067)

Lane 2: IP with Anti-p62 C-terminal pAb (MBL; code no. PM066)

Immunoblotted with Anti-p62 (SQSTM1) pAb (MBL; code no. PM045)

**Flow cytometric analysis for adherent cells**

- 1) Detach the cells from culture dish.
- 2) Wash the cells 3 times with 1 mL of washing buffer [PBS containing 2% fetal calf serum (FCS)].
- 3) Add 200 µL of 4% paraformaldehyde (PFA) to the cell pellet after tapping. Mix well, then fix the cells for 10 min. at room temperature.
- 4) Wash the cells 2 times with 1 mL of washing buffer.
- 5) Add 200 µL of PBS containing 100 µg/mL Digitonin to the cell pellet after tapping. Mix well, then permeabilize the cells for 10 min. at room temperature.
- 6) Wash the cells 2 times with 1 mL of washing buffer.
- 7) Resuspend the cells with washing buffer ( $5 \times 10^6$  cells/mL).
- 8) Add 100 µL of the cell suspension into each tube, and centrifuge at 500 x g for 1 min. at room temperature (20~25°C). Remove the supernatant by careful aspiration.
- 9) Add 20 µL of Clear Back (human Fc receptor blocking reagent, MBL; code no. MTG-001) to the cell pellet after tapping. Mix well and incubate for 5 min. at room temperature.
- 10) Add the isotype control antibody at the concentrations comparable to those of the specific antibody of interest. Mix well and incubate for 30 min. at room temperature.
- 11) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 min. at room temperature. Remove supernatant by careful aspiration.
- 12) Add FITC-conjugated anti-Guinea Pig IgG antibody diluted with the washing buffer. Mix well and incubate for 30 min. at room temperature.
- 13) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 min. at room temperature. Remove supernatant by careful aspiration.
- 14) Resuspend the cells with 500 µL of the washing buffer and analyze by a flow cytometer.



***Flow cytometric analysis of p62 in HeLa***

Closed: Secondary antibody alone

Open: Normal Guinea Pig IgG (PM067) or Anti-p62 C-terminal pAb (MBL; code no. PM066)