## T-Select

# Human CD1d Tetramer

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

## Background :

Natural Killer T (NKT) cells, a type of T cell that plays a significant role in immune response, produces a large quantity of INF- $\gamma$  and IL-4 in response to glycolipids that presented by CD1d molecules. In recent years, NKT cells are reported to play a part in diabetes and tumor immunity. Therefore, a technology allowing quantitative measurement of CD1d-positive NKT cells would be a useful tool for immunology and clinical laboratory examinations.

The development of MHC Tetramer technology has provided a breakthrough in the ability to follow T cell populations defined by their antigen specificity. Tetramers have been used widely to obtain a detailed analysis of the distribution and frequency of conventional CD4<sup>+</sup> and CD8<sup>+</sup> antigen-specific T cells during a variety of immune responses. T-Select Human CD1d Tetramer is a reagent created by tetramerizing biotinylated human CD1d/ $\beta$ 2m complexes with streptavidin labeled phycobiliprotein.  $\alpha$ -galactosylceramide ( $\alpha$ -GalCer), a glycosphingolipid originally isolated from marine sponges, appears to be presented by CD1d to activate both human and mouse NKT cells.  $\alpha$ -GalCer loaded T-Select Human CD1d Tetramer is a highly specific reagent for detection of NKT cells. Measurement can be performed using isolated lymphocytes/monocytes.

## Specificity:

T-Select Human CD1d Tetramer recognizes human NKT cells that bind specifically with CD1d molecule and  $\alpha$ -GalCer complex.

## Storage:

Store protected from light at 2-8°C. Do not freeze.

## Fluorescent labeling :

Phycoerythrin

Excitation wavelength: 486 - 580 nm

Fluorescent wavelength: 586 - 590 nm

Allophycocyanin

Excitation wavelength: 633 - 635 nm

Fluorescent wavelength: 660 - 680 nm

## Formulation:

The reagent is dissolved in aqueous buffer that contains 10 mM Tris-HCl (pH 8.0), 150 mM NaCl, 0.5 mM EDTA, 0.2% BSA and 0.09% NaN<sub>3</sub>.

\*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.

## Stability:

The product is stable for at least 2 year when stored as recommended. If the reagent shows physical change (usually it is a slightly pinkish or blue liquid), do not use the reagent because it may be degraded.

## Preparation of reagent:

T-Select Human CD1d Tetramer does not contain  $\alpha$ -GalCer. Please combine with  $\alpha$ -GalCer before use according to the procedure described below.

1. Dissolve  $\alpha$ -GalCer in pyridine to make a concentration of 1 mg/mL. Use this as a stock solution, and dilute it with 0.5% Tween-20, 0.9% NaCl to make a concentration of 200  $\mu$ g/mL.
2. Add 5  $\mu$ L of diluted  $\alpha$ -GalCer to 100  $\mu$ L of T-Select Human CD1d Tetramer-PE, and pipette the whole solution slowly. Incubate it at room temperature (20-30°C) or 37°C for 12 to 18 hours. Minimize exposure to light.
3. After that, store the reagent as is at 4°C.

## Cell staining procedure:

1. Prepare single cell suspension from anticoagulated human peripheral blood according to the standard procedure. For staining, suspend the cells in FCM buffer (2% FCS/0.05% NaN<sub>3</sub>/PBS) at a concentration of up to  $1 \times 10^7$  cells/mL.
2. Add 100  $\mu$ L of cell suspension to each test tube and centrifuge it at 400 x g for 5 minutes. Aspirate the supernatant.
3. Add 10  $\mu$ L of Clear Back (human Fc receptor blocking reagent, MBL code no. MTG-001) and 30  $\mu$ L of FCM buffer to the cell pellet after tapping. Mix well and incubate for 5 minutes at room temperature.
4. Add 10  $\mu$ L of T-Select Human CD1d Tetramer-SA-PE to each test tube and mix well. Incubate the cells in the dark for 20 minutes at 4°C.
5. Add 1 mL of FCM buffer.
6. Centrifuge at 400 x g for 5 minutes. Aspirate the supernatant. Repeat the same operation for three times.
7. Suspend the pellet with 500  $\mu$ L of FCM buffer. Analyze it immediately, or suspend it with 0.5% paraformaldehyde/PBS, and store the sample in dark room at 4°C. Be sure to analyze it within 24 hours.

\* To add antibodies such as CD3, add them at the step 4.



## Consideration :

- A. If the blood erythrocyte remains in the cell sample, we recommend hemolyzing them. If the blood erythrocyte still remains after being hemolyzed, we recommend staining the cells with anti-CD45 antibody simultaneously and analyzing the result with lymphocyte gating.
- B. We recommend using Clear Back (MBL code no. MTG-001) to reduce nonspecific staining of cells by endocytosis in the macrophages.
- C. For staining of *in vitro* cultured lymphocyte, we recommend staining with 7-AAD for exclusion of dead cells and non-viable cells.
- D. Paraformaldehyde fixation of cells is not needed if the cells are analyzed within a couple of hours after staining.

**Precautions:**

1. The reagent contains 0.09% of sodium azide. Sodium azide under acid conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be flushed with running water while being discarded. These precautions are recommended to avoid deposits in metal piping in which explosive conditions can develop. If skin or eye contact occurs, wash excessively with water.
2. Samples and all materials coming in contact with them should be handled as if capable of transmitting infection and disposed of with proper precautions.
3. Never pipet by mouth and avoid contact of samples with skin and mucous membranes.
4. Minimize exposure of reagents to light during storage or incubation.
5. Recommended cell survival rate of venous blood specimen is  $\geq 90\%$ .
6. To obtain appropriate result with whole blood, we recommend to keeping the test sample in a blood collection tube at room temperature and turning upside-down repeatedly just before staining. Do not use cold test blood in order to have appropriate result.
7. Do not incubate the cells for a long time with the hemolysis reagent. The long incubation results in the disruption of leukocyte.
8. Erythrocyte of abnormal test blood, such as nucleated erythrocyte, blood of abnormal hemoglobin disease, can not be hemolyzed well. In such case, the unhemolyzed erythrocyte is improperly counted as a leukocyte. This improper count results in increasing the number of leukocyte and decreasing the number of the positive rate of NKT cells.

Example of human PBMC staining using T-Select Human CD1d Tetramer-PE with (upper) or without  $\alpha$ -GalCer (lower)

**References:**

- 1) Matsuda JL, *et al. J. Exp. Med.* **192**: 741-754 (2000)
- 2) Karadimitris A, *et al. PNAS* **98**: 3294-3298 (2001)
- 3) Kita H, *et al. Gastroenterology* **123**: 1031-1043 (2002)
- 4) Sidobre S, and Kronenberg M, *J. Immunol. Methods* **268**: 107-121 (2002)
- 5) Wu D, *et al. PNAS* **102**: 1351-1356 (2005)
- 6) Li D, *et al. J. Immunol.* **182**: 1033-1040 (2009)

**Related Products:**

**CD1d Tetramers**

TS-HCD-1	Human CD1d Tetramer-PE (50 tests)
TS-HCD-2	Human CD1d Tetramer-APC (50 tests)
TS-MCD-1	Mouse CD1d Tetramer-PE (50 tests)
TS-MCD-1S	Mouse CD1d Tetramer-PE (10 tests)
TS-MCD-2	Mouse CD1d Tetramer-APC (50 tests)

**Antibody for NKT**

IM-1588	anti-TCR V $\alpha$ 24	(clone C15)
IM-1589	FITC labeled anti-TCR V $\alpha$ 24	(clone C15)
IM-2283	PE labeled anti-TCR V $\alpha$ 24	(clone C15)
IM-2027	Biotin labeled anti-TCR V $\alpha$ 24	(clone C15)
IM-1582	anti-TCR V $\beta$ 11	(clone C21)
IM-1586	FITC labeled anti-TCR V $\beta$ 11	(clone C21)
IM-2290	PE labeled anti-TCR V $\beta$ 11	(clone C21)
IM-2018	Biotin labeled anti-TCR V $\beta$ 11	(clone C21)

6604629	CD3	(clone UCHT1)
A07746	CD3-FITC	(clone UCHT1)
A07747	CD3-PE	(clone UCHT1)
A07749	CD3-PC5	(clone UCHT1)
IM-2467	CD3-APC	(clone UCHT1)

IM-0398	CD4	(clone 13B8.2)
A07750	CD4-FITC	(clone 13B8.2)
A07751	CD4-PE	(clone 13B8.2)
A07752	CD4-PC5	(clone 13B8.2)
IM-2468	CD4-APC	(clone 13B8.2)

6602139	CD8	(clone T8)
6603861	CD8-FITC	(clone T8)
A07757	CD8-PE	(clone B9.11)
6607011	CD8-PC5	(clone T8)
IM-2469	CD8-APC	(clone B9.11)

A07782	CD45-FITC	(clone J.33)
A07783	CD45-PE	(clone J.33)
A07785	CD45-PC5	(clone J.33)
IM-2473	CD45-APC	(clone J.33)

IM-0780	CD69	(clone TP1.55.3)
IM-1943	CD69-PE	(clone TP1.55.3)
IM-2656	CD69-PC5	(clone TP1.55.3)

IM-1610	CD94	(clone HP-3B1)
IM-2276	CD94-PE	(clone HP-3B1)

K0061-3	CD161	(clone HP-3G10)
K0061-4	CD161-FITC	(clone HP-3G10)

**Antibody for NK**

6604894	CD16-FITC	(clone 3G8)
A07766	CD16-PE	(clone 3G8)
A07767	CD16-PC5	(clone 3G8)
6607118	CD16-PC7	(clone 3G8)

6602705	CD56	(clone N901)
A07788	CD56-PE	(clone N901)
A07789	CD56-PC5	(clone N901)
A21692	CD56-PC7	(clone N901)
IM-2474	CD56-APC	(clone N901)

IM-1166	CD57	(clone NC1)
IM-0466	CD57-FITC	(clone NC1)

**Others**

MTG-001	Clear Back (Human FcR blocking reagent)
A07704	7-AAD Viability Dye
IM-1981	IL-4 EIA Kit
IM-1743	IFN- $\gamma$ EIA Kit
TS-8005	T-Select MHC IFN- $\gamma$ Kit
IM-1400	OptiLyse B
A11895	OptiLyse C

Please check our web site (<http://ruo.mbl.co.jp>) for up-to-date information on products and custom MHC Tetramers.