



Honenpak C18 (for ABEE Labeling Kit)

Pkg.	1 piece (75 mm X 4.6mm I.D.)
Usable pH range	pH 1-9
Feature	This HPLC column has been developed specifically for the analysis of sugars labeled with the ABEE Labeling Kit (1,2). The column is shipped filled with a 70% acetonitrile solution.
Storage	After replacement with a 50% or higher acetonitrile solution, the column should be sealed tightly and stored in a place with low temperature variation.

Reference

- 1) Yasuno, S., et al., *Biosci. Biotechnol. Biochem.*, 61, 1944 (1997)
- 2) Yasuno, S., et al., *Biosci. Biotechnol. Biochem.*, 63, 1353 (1999)

For reserach use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioua.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioua.com