

FITC-labeled Gelatin Zymogram Precast Gel, Cosmo Bio-type

PMC-AK84-COS

PMC-AK85-COS

Last Updated on Oct. 1st, 2016

www.cosmobio.com

These are Precast Gel for FITC-labeled Gelatin Zymography Kit (PMC-AK83-COS).

【 1 】 Component

Cat No.	Description	Quantity	Storage
PMC-AK84-COS	Precast gel for 12 samples	5 pieces	4°C, Do not freeze
PMC-AK85-COS		10 pieces	

The size of the glass plate gel is 100 mm (W) × 100 mm (H) and the thickness is 1mm.

【 1-2 】 Features of FITC-labeled Gelatin-zymography Kit

MMPs are widely detected by zymography. The FITC-labeled Gelatin-Zymography Kit (Cat.No.AK83-COS) provides an easy system of the electrophoresis for zymography. This product is used for detecting ProMMP-2, MMP-2 and ProMMP-9 in blood, body fluid, secretion, cell lysate, cell culture medium, and other samples.

Cat No.	Description	Components
PMC-AK83-COS	FITC-labeled Gelatin-zymography Kit (Cosmo Bio Type)	5 pieces of Precast gel for 12 samples 2 bottles of 100 mL Electrophoresis Buffer (× 10) 100 mL Washing Buffer (× 10) 25 mL Reaction Buffer (× 10) 5 mL Sample Preparation Buffer 0.2 mL MMP markers (ProMMP-2, MMP-2, ProMMP-9)
PMC-AK38-COS	MMP markers	200 μL ((ProMMP-2, MMP-2, ProMMP-9)



COSMO BIO Co., LTD.

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
TEL : (81)3-5632-9617
FAX : (81)3-5632-9618
e-mail : export@cosmobio.co.jp
URL : www.cosmobio.com



www.cosmobio.com