

T-Select MHC Tetramer

# HLA-C\*03:03 EBV EBNA1

## Tetramer-FVYGGSKTSL (50 tests)

For Research Use Only. Not for use in diagnostic procedures.  
MBL manufactures and distributes these products under license from Beckman Coulter, Inc..

### Background

T lymphocytes play a central role in immune system. Total T cell and T cell subset counts are measured by detection of various cell surface molecules. Enumeration of CD8<sup>+</sup> antigen-specific T cells requires cognate recognition of the T cell receptor (TCR) by a class I MHC/peptide complex. This can be done using T-Select MHC class I Tetramers which are composed of four MHC class I molecules each bound to the specific peptide and conjugated with a fluorescent protein. Thus, T-Select MHC Tetramer assays allow quantitation of the total T cell population specific for a given peptide complexed in a particular MHC molecule. Furthermore, since binding does not depend on functional pathways, this population includes specific CD8<sup>+</sup> T cells regardless of functional status. Measurements may be performed in whole blood or isolated lymphocyte/mononuclear cell preparations. In some cases where frequency is low, it may be necessary to perform an *in vitro* cell expansion. Specific cell staining is accomplished by incubating the sample with the T-Select MHC Tetramer reagent, then washing away excess Tetramer. The number of Tetramer positive lymphocytes is then determined by flow cytometry.

Epstein-Barr virus (EBV) is a latent  $\gamma$ -herpesvirus with high prevalence in human population worldwide. EBV was the first human virus found to be related with tumor. It is closely associated with a number of tumors such as Burkitt's lymphoma, Hodgkin's lymphoma, nasopharyngeal carcinoma (NPC) and some gastric carcinoma. Furthermore, recent studies suggested that EBV might also contribute to certain autoimmune diseases.

Because of its oncogenic importance, EBV targeted immunotherapies are widely studied. The best characterized case is EBV-induced post-transplant lymphoproliferative disorders (PTLD), which continues to be a major complication after solid organ or bone marrow transplantation in high-risk patients.

Several studies have shown that adoptive transfer of donor-derived EBV-specific T cells is an effective strategy for both prevention and treatment of PTLD. For other EBV-positive tumors, such as NPC, autologous CTL therapy has been evaluated in many clinical trials, with promising results.

These Tetramer reagents comprise human class I HLA-C\*03:03 and epitope peptide derived from Epstein-Barr nuclear antigen 1 (EBNA1), and they can detect HLA-C\*03:03-restricted EBV EBNA1-specific CD8<sup>+</sup> T cells by flow cytometry.

EBNA1 is the only latent viral protein known to be consistently expressed across all EBV-associated tumors. Therefore, EBNA1 has long been considered an ideal target for immunotherapy.

### HLA Restriction

HLA-C\*03:03

### Origin and Sequence of CTL Epitope

EBV EBNA1 (508-517 aa, FVYGGSKTSL)

### Reagents

500  $\mu$ L liquid - 10  $\mu$ L/test

The Tetramer is dissolved in an aqueous buffer containing 0.5 mM EDTA, 0.2% BSA, 10 mM Tris-HCl (pH 8.0), 150 mM NaCl, and 0.09% NaN<sub>3</sub>.

### Conjugates

TS-M150-1

Streptavidin-Phycoerythrin (SA-PE)

Excites at 486-580 nm

Emits at 586-590 nm

TS-M150-2

Streptavidin-Allophycocyanin (SA-APC)

Excites at 633-635 nm

Emits at 660-680 nm

### Storage Conditions

Store at 2 to 8°C. Do not freeze. Minimize exposure to light. The expiration date is indicated on the vial label.

### Evidence of Deterioration

Any change in the physical appearance of this reagent may indicate deterioration and the reagent should not be used. The normal appearance is a clear, colorless to pink (SA-PE), or light blue (SA-APC).

## Usage

This reagent is for use with standard flow cytometry methodologies.

## References for This Product

- 1) Ito Y, *et al. J Gen Virol* **88**: 770–780 (2007)
- 2) Marescotti D, *et al. Immunology* **129**: 386-395 (2010)
- 3) Chang CX, *et al. Eur J Immunol* **43**: 1109-1120 (2013)
- 4) Westhoff SD, *et al. Viruses* **5**: 226-240 (2013)

## High Specificity

The T cell surface CD8 enhances T cell antigen recognition by binding to HLA class I molecules. Therefore, MBL produced T-Select MHC class I human Tetramers with one point mutation at the HLA  $\alpha$ 3 domain known to alter the interaction with CD8. These mutated Tetramers showed a greatly diminished nonspecific binding but retained specific binding. Alterations of CD8 binding by mutation of the MHC greatly improved the specificity of MHC-peptide multimers, thus providing efficient tools to sort specific human T cells for immunotherapy. (French application Number; FR9911133)

## References for T-Select MHC Tetramer

- Altman JD, *et al. Science* **274**: 94-96 (1996)  
McMichael AJ, *et al. J Exp Med* **187**: 1367-1371 (1998)  
Bodinier M, *et al. Nat Med* **6**: 707-710 (2000)

## Statement of Warnings

1. This reagent contains 0.09% sodium azide. Sodium azide under acid conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be flushed with running water while being discarded. These precautions are recommended to avoid deposits in metal piping in which explosive conditions can develop. If skin or eye contact occurs, wash excessively with water.
2. Specimens, samples and material coming in contact with them should be handled as if capable of transmitting infection and disposed of with proper precautions.
3. Never pipette by mouth and avoid contact of samples with skin and mucous membranes.
4. Minimize exposure of reagent to light during storage or incubation.
5. Avoid microbial contamination of reagent or erroneous results may occur.
6. Use Good Laboratory Practices (GLP) when handling this reagent.

## Materials Required But Not Supplied

- 12 x 75 mm polypropylene test tubes
- Transfer pipettes
- Pipettors and disposable pipette tips
- Vortex mixer
- Centrifuge capable of 150 x g or 400 x g
- Aspirator

- PBS
- Red blood cell lysis reagent
- Anti-CD8-FITC, Beckman Coulter, Inc., PN 6603861
- Anti-CD8-PC5, Beckman Coulter, Inc., PN 6607011
- 7-AAD Viability Dye, Beckman Coulter, Inc., PN A07704
- Clear Back (human FcR blocking reagent), MBL, PN MTG-001

## Procedure for Whole Blood

1. Collect blood by venipuncture into a blood collection tube containing an appropriate anti-coagulant.
2. Add 10  $\mu$ L of T-Select MHC Tetramer to each 12 x 75 mm test tube.
3. Add 200  $\mu$ L of whole blood into each test tube.
4. Vortex gently.
5. Incubate for 30-60 minutes at 2-8°C or room temperature (15-25°C) protected from light.
6. Add any additional antibodies (e.g. anti-CD8) and vortex gently.
7. Incubate for 30 minutes at 2-8°C protected from light.
8. Lyse red blood cells using commercially available reagents.
9. Prepare samples according to description of the package insert.
10. Store prepared samples at 2-8°C protected from light for a minimum of 1 hour (maximum 24 hours) prior to analysis by flow cytometry.

## Procedure for Peripheral Blood Mononuclear Cells

1. Prepare peripheral blood mononuclear cells (PBMC) according to established procedures. Cells should be re-suspended at a concentration of  $2 \times 10^7$  cells/mL. 50  $\mu$ L of sample is required for each T-Select MHC Tetramer determination.
2. Add 10  $\mu$ L of Clear Back (human FcR blocking reagent, MBL, PN MTG-001) to each 12 x 75 mm test tube.
3. Add 50  $\mu$ L PBMC into each test tube (e.g.  $1 \times 10^6$  cells per tube).
4. Incubate for 5 minutes at room temperature.
5. Add 10  $\mu$ L of T-Select MHC Tetramer and vortex gently.
6. Incubate for 30-60 minutes at 2-8°C or room temperature (15-25°C) protected from light.
7. Add any additional antibodies (e.g. anti-CD8) and vortex gently.
8. Incubate for 30 minutes at 2-8°C protected from light.
9. Add 3 mL of PBS or FCM buffer (2% FCS/0.09%  $\text{NaN}_3$ /PBS).
10. Centrifuge tubes at 400 x g for 5 minutes.
11. Aspirate or decant the supernatant.
12. Resuspend the pellet in 500  $\mu$ L of PBS with 0.5% formaldehyde.
13. Analyze prepared samples by flow cytometry. If necessary, store the samples at 2-8°C protected from light for a maximum of 24 hours prior to analysis.

## Limitations

1. For optimal results with whole blood, retain specimens in blood collection tubes at room temperature, while rocking, prior to staining and analyzing. Refrigerated specimens may give aberrant results.
2. Recommended cell viability for venous blood specimens is > 90%.
3. Prolonged exposure of cells to lytic reagents may cause white blood cell destruction and loss of cells in the population of interest.
4. All red blood cells may not lyse under the following conditions: nucleated red blood cells, abnormal protein concentration or hemoglobinopathies. This may cause falsely decreased results due to unlysed red blood cells being counted as leukocytes.

## Technical Hints

- A. If PBMC culture is needed, we recommend the use of heparin as an anti-coagulant.
- B. Clear Back reagent (human FcR blocking reagent) may effectively block non-specific binding caused by macrophages or endocytosis, resulting in clear staining when cells are stained with MHC Tetramer and antibodies. Please refer to the data sheet (MBL, PN MTG-001) for details.
- C. A Tetramer that is constructed with the same allele of interest and an irrelevant peptide may be used as a negative control.
- D. We recommend the use of CD8 antibody, clone SFC121Thy2D3 (T8, Beckman Coulter, Inc.), which does not block or interfere with the specific binding of MHC Tetramers to T cells.
- E. The use of CD45 antibody and gating of the lymphocyte population are recommended in order to reduce contamination of unlysed or nucleated red blood cells in the gate.
- F. Apoptotic, necrotic, and/or damaged cells are sources of interference in the analysis of viable cells by flow cytometry. Cell viability should be determined by 7-aminoactinomycin D (7-AAD) staining; intact viable cells remain unstained (negative).
- G. Cells do not require fixation prior to analysis if the stained cells are analyzed by flow cytometry within several hours.

## Related Products

### T-Select Human Tetramers

#### EBV

TS-0011-1C	HLA-A*02:01 EBV BMLF1 Tetramer-GLCTLVAML-PE
TS-M006-1	HLA-A*02:01 EBV LMP1 Tetramer-YLQQNWWTL-PE
TS-M030-1	HLA-A*02:01 EBV LMP2 Tetramer-TVCGGIMFL-PE
TS-M031-1	HLA-A*02:01 EBV LMP2 Tetramer-LLWTLVLL-PE
TS-M069-1	HLA-A*02:01 EBV LMP2 Tetramer-FLYALALL-PE
TS-M032-1	HLA-A*02:01 EBV LMP2 Tetramer-CLGGLLTMV-PE
TS-M003-1	HLA-A*24:02 EBV BMLF1 Tetramer-DYNFVKQLF-PE
TS-M002-1	HLA-A*24:02 EBV BRLF1 Tetramer-TYPVLEEMF-PE

TS-M004-1	HLA-A*24:02 EBV EBNA3A Tetramer-RYSIFFDYM-PE
TS-M005-1	HLA-A*24:02 EBV EBNA3B Tetramer-TYSAGIVQI-PE
TS-M034-1	HLA-A*24:02 EBV LMP2 Tetramer-PYLFWLAAI-PE
TS-M001-1	HLA-A*24:02 EBV LMP2 Tetramer-IYVLVMLVL-PE
TS-M035-1	HLA-A*24:02 EBV LMP2 Tetramer-TYGPVFMSL-PE
TS-M009-1	HLA-A*24:02 EBV Mix Tetramer-PE
TS-M028-1	HLA-A*11:01 EBV EBNA3B <sub>399-408</sub> Tetramer-PE
TS-M029-1	HLA-A*11:01 EBV EBNA3B <sub>416-424</sub> Tetramer-PE
TS-M111-1	HLA-A*11:01 EBV LMP2 Tetramer-SSCSCPLSK-PE
TS-M135-1	HLA-A*11:01 EBV LMP2 S9T Tetramer-SSCSCPLTK-PE
TS-M036-1	HLA-A*08:01 EBV BZLF1 Tetramer-RAKFKQLL-PE
TS-M037-1	HLA-B*35:01 EBV BZLF1 Tetramer-EPLPQQQLTAY-PE
TS-M038-1	HLA-B*35:01 EBV LMP2 Tetramer-MGSLEMVPM-PE
TS-M142-1	HLA-B*07:02 EBV EBNA3A Tetramer-RPPIFIRRL-PE
TS-M803-1	HLA-DRB1*01:01 EBV EBNA1 <sub>515-527</sub> Tetramer-PE

## Control

TS-M007-1	HLA-A*24:02 Negative Tetramer-RYLDRDQQLL-PE
TS-M007-2	HLA-A*24:02 Negative Tetramer-RYLDRDQQLL-APC
TS-M007-3	HLA-A*24:02 Negative Tetramer-RYLDRDQQLL-FITC
TS-0029-1C	HLA-A*02:01 Negative Tetramer-PE
TS-0029-2C	HLA-A*02:01 Negative Tetramer-APC
TS-M801-1	HLA-DRB1*01:01 human CLIP <sub>103-117</sub> Tetramer-PE
TS-M802-1	HLA-DRB1*01:01 HIV gag <sub>295-307</sub> Tetramer-PE

## T-Select PEPTIDES

TS-M001-P	HLA-A*24:02 EBV LMP2 peptide
TS-M002-P	HLA-A*24:02 EBV BRLF1 peptide
TS-M003-P	HLA-A*24:02 EBV BMLF1 peptide
TS-M004-P	HLA-A*24:02 EBV EBNA3A peptide
TS-M005-P	HLA-A*24:02 EBV EBNA3B peptide
TS-M035-P	HLA-A*24:02 EBV LMP2 <sub>419-427</sub> peptide
TS-M007-P	HLA-A*24:02 HIV env gp160 peptide
TS-0011-P	HLA-A*02:01 EBV BMLF1 peptide
TS-M069-P	HLA-A*02:01 EBV LMP2 <sub>356-364</sub> peptide
TS-0029-P	HLA-A*02:01 Negative peptide
TS-M111-P	HLA-A*11:01 EBV LMP2 peptide
TS-M135-P	HLA-A*11:01 EBV LMP2 S9T peptide
TS-M803-P	HLA-DRB1*01:01 EBV EBNA1 <sub>515-527</sub> Peptide
TS-M801-P	Human CLIP <sub>103-117</sub> peptide
TS-M802-P	HLA-DRB1*01:01 HIV gag <sub>295-307</sub> peptide

## Kits

4844	IMMUNOCYTO CD107a Detection Kit
AM-1005M	IMMUNOCYTO Cytotoxicity Detection Kit

## Others

6603861	CD8-FITC (T8)
6607011	CD8-PC5 (T8)
A07704	7-AAD Viability Dye
MTG-001	Clear Back (Human FcR blocking reagent)

Please check our web site (<http://ruo.mbl.co.jp>) for up-to-date information on products and custom MHC Tetramers.

### Experimental Data

The frozen stock of PBMCs of healthy HLA-C\*03:03-positive donors were thawed and evaluated.

Approximately  $3 \times 10^5$  cells were stained with MHC Tetramers, CD8 antibody, and 7-AAD.

The lymphocyte population was defined by an FSC/SSC gate, and the viable cell population was defined by an FSC/7-AAD. Data were analyzed by double gating on the lymphocyte and viable cell population. Numbers in the top right quadrants represent the percentage of MHC Tetramer-positive cells in the total CD8<sup>+</sup> cells.

As a result, HLA-C\*03:03-restricted EBNA1-specific CTLs were detected in 2 donors.

