



RANDOX *life sciences*

PRODUCT SPECIFICATION

Product	:	Anti- Cyclin-Dependent Kinase 4 (CDK4)
Cat. No.	:	PAS9321
Lot No.	:	RS415-17
Immunogen	:	Recombinant Human CDK4
Host Animal	:	Sheep
Format	:	Affinity purified Ig
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	≤ -20°C

Randox Life Sciences

55 Diamond Road, Crumlin, Co. Antrim, United Kingdom, BT29 4QY
t: +44 (0) 289442 2413 f: +44 (0) 289445 2912 e: lifesciences@randox.com
www.randox-lifesciences.com