



PRODUCT SPECIFICATION

Product	:	Anti-Chloramphenicol
Cat. No.	:	PAS 9511
Lot No.	:	RS485-15
Immunogen	:	Chloramphenicol-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	9.88mg/ml (U.V. abs @ 280nm)
Working Concentration	:	10µg/ml (EIA)
pH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

Compound	% Cross Reactivity
Chloramphenicol	100
Thiamphenicol	4
Florphenicol	15

Randox Life Sciences