



PRODUCT SPECIFICATION

Product	:	Anti-17 α -Trenbolone
Cat. No.	:	PAS 9604
Lot No.	:	RS729-08
Immunogen	:	17 α -Trenbolone(3)-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	9.59mg/ml (U.V. abs @ 280nm)
Working Concentration	:	0.5 μ g/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

Compound	Cross Reactivity
Trenbolone (17 β)	1.41%
Trenbolone (17 α)	266.14%
Methyltestosterone	<1%
Boldenone (17 α)	1.80%
Boldenone (17 β)	<1%
Methandriol	<1%
CLAD	<1%
Medroxyprogesterone acetate	<1%
Ethinyl Estradiol	<1%
Stanozolol	<1%

Randox Life Sciences