



PRODUCT SPECIFICATION

Product	:	Anti-Clenbuterol
Cat. No.	:	PAS9683
Lot No.	:	104-03
Immunogen	:	Clenbuterol(N)-HSA
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	10mg/ml (U.V. abs @ 280nm)
Working Concentration	:	10µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

SPECIFICITY:

Compound	% Cross-Reactivity
Clenbuterol	100
Mabuterol	114
Brombuterol	96
Methyl-Clenbuterol	53
Mapenterol	49
Salbutamol	3.8
Terbutaline	1.6
Cimaterol	< 1.0
Pirbuterol	0.9
Ractopamine	< 0.2
Feneterol	< 0.2

Randox Life Sciences