



RANDOX *life sciences*

PRODUCT SPECIFICATION

| | | |
|------------------------------|---|---------------------------------------|
| Product | : | Anti-Creatinine |
| Cat. No. | : | PAS 9686 |
| Lot No. | : | RS300-24 |
| Immunogen | : | Creatinine(N)-BSA |
| Host Animal | : | Sheep |
| Format | : | Ig fraction |
| Ig Concentration | : | 6.13mg/ml (U.V. abs @ 280nm) |
| Working Concentration | : | 2µg/ml (EIA) |
| PH | : | 7.2 |
| Buffer | : | 20mM Phosphate, 150mM Sodium Chloride |
| Preservative | : | 0.09% Sodium Azide |
| Storage | : | -20°C |

Randox Life Sciences

55 Diamond Road, Crumlin, Co. Antrim, United Kingdom, BT29 4QY
t: +44 (0) 289442 2413 f: +44 (0) 289445 2912 e: lifesciences@randox.com
www.randox-lifesciences.com