



RANDOX *life sciences*

PRODUCT SPECIFICATION

Product	:	Anti-Cyclic Adenosine Monophosphate (Cyclic AMP)
Cat. No.	:	PAS 9687
Lot No.	:	RS150-16
Immunogen	:	Cyclic AMP-BSA
Host Animal	:	Sheep
Format	:	Ig fraction
Ig Concentration	:	6.05mg/ml (U.V. abs @ 280nm)
Working Concentration	:	6.05µg/ml (EIA)
pH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	-20°C

Randox Life Sciences

55 Diamond Road, Crumlin, Co. Antrim, United Kingdom, BT29 4QY
t: +44 (0) 289442 2413 f: +44 (0) 289445 2912 e: lifesciences@randox.com
www.randox-lifesciences.com