



RANDOX *life sciences*

PRODUCT SPECIFICATION

Product	:	Anti-Androstenedione
Cat. No.	:	PAS9853
Lot No.	:	RS1297-B3
Immunogen	:	Androstenedione(3)-BTG
Host Animal	:	Sheep
Format	:	Ig fraction
Ig concentration	:	7.03mg/ml (U.V.abs @ 280nm)
Working Concentration	:	0.625µg/ml (EIA)
PH	:	7.2
Buffer	:	20mM Phosphate, 150mM Sodium Chloride
Preservative	:	0.09% Sodium Azide
Storage	:	≤-20°C

Randox Life Sciences

55 Diamond Road, Crumlin, Co. Antrim, United Kingdom, BT29 4QY
t: +44 (0) 289442 2413 f: +44 (0) 289445 2912 e: lifesciences@randox.com
www.randox-lifesciences.com