



Anti Cathepsin L2 [Clone : 236-8C9]

Product type	Monoclonal antibody
Immunogen	Recombinant human Cathepsin L2.
Host	Mouse
Clone number	236-8C9
Isotype	IgG2b/k
Purification	Purified
Buffer	100 µg purified antibody lyophilized from 0.1 M Na-Phosphate buffer, pH7.0 with 1 mg of protease-free BSA. We recommend resuspending the lyophilized antibody with 0.2 mL distilled water.
Volume	100 µg
Label	Unlabeled
Specificity	The human Cathepsin L2, but not with human Cathepsins L, Z and F on western blotting.
Reactivity	Human. Other species are not tested.
Storage	Store F-106 at 4 centigrade until reconstituted, then store in aliquots at -20 centigrade or at 4 centigrade with 0.05% azide. Freezing of aliquots is best for long term storage of reconstituted product; repetitive freezing and thawing should be avoided.

Application notes

Recommended dilutions
• Western Blotting : 2.0 µg/mL

Other applications have not been tested.
Optimal dilutions/concentrations should be determined by the end user.

References 1) Santamaria, I. et al., Cancer Res., 58, 1624-1630, 1998.



COSMO BIO Co., LTD.

【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com