



## Anti Cathepsin Z [Clone : 237-11B4]

<b>Product type</b>	Monoclonal antibody
<b>Immunogen</b>	Recombinant human Cathepsin Z
<b>Host</b>	Mouse
<b>Clone number</b>	237-11B4
<b>Isotype</b>	IgG1/k
<b>Purification</b>	Purified
<b>Buffer</b>	100 µg purified antibody lyophilized from 0.1 M Na-Phosphate buffer, pH7.0 with 1 mg of protease-free BSA. We recommend resuspending the lyophilized antibody with 0.2 mL distilled water.
<b>Volume</b>	100 µg
<b>Label</b>	Unlabeled
<b>Specificity</b>	The human Cathepsin Z, but not with human Cathepsins L, L2 and F on western blotting.
<b>Reactivity</b>	Human. Other species are not tested.
<b>Storage</b>	Store F-108 at 4 centigrade until reconstituted, then store in aliquots at -20 centigrade or at 4 centigrade with 0.05% azide. Freezing of aliquots is best for long term storage of reconstituted product; repetitive freezing and thawing should be avoided.

### Application notes

Recommended dilutions  
• Western Blotting : 2.0 µg/mL

Other applications have not been tested.  
Optimal dilutions/concentrations should be determined by the end user.

**References** 1) Santamaria, I. et al., J. Biol. Chem., 273, 16816-16823, 1998.



COSMO BIO Co., LTD.

【JAPAN】  
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】  
2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobiousa.com](mailto:info@cosmobiousa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobiousa.com](http://www.cosmobiousa.com)