



Anti TIMP-4 [Clone : 216-2C9]

Product type	Monoclonal antibody
Immunogen	Recombinant human tissue inhibitor of metalloproteinases-4.
Host	Mouse
Clone number	216-2C9
Isotype	IgG2b/k
Purification	Purified
Buffer	0.1 M Na-Phosphate buffer, pH7.0 containing 0.5% protease-free bovine serum albumin and 0.05% sodium azide.
Concentration	500 µg/mL
Volume	200 µL
Label	Unlabeled
Specificity	hTIMP-4 and does not cross-react with hTIMP-1, hTIMP-2 and hTIMP-3.
Reactivity	Human. Other species are not tested.
Storage	Maintain at 2-8°C in undiluted aliquots for up to 12 months. Do not freeze.

Application notes

Recommended dilutions

- Western Blotting : 5 µg/mL

Other applications have not been tested.

Optimal dilutions/concentrations should be determined by the end user.

References

- 1) Greene, J. et al., J. Biol. Chem., 271, 30375-30380, 1996.
- 2) Douglas, D. A. et al., J. Protein Chem., 16, 237-255, 1997.
- 3) Liu, Y. E. et al., J. Biol. Chem., 272, 20479-20483, 1997.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com