



Anti Urokinase (H) [Clone : 5-15F12]

| | |
|----------------------|---|
| Product type | Monoclonal antibody |
| Immunogen | Human Urokinase |
| Host | Mouse |
| Clone number | 5-15F12 |
| Isotype | IgG1/k |
| Purification | Purified |
| Buffer | 0.1 M Na-Phosphate buffer, pH7.0 containing 0.5% protease-free bovine serum albumin and 0.05% sodium azide. |
| Concentration | 500 µg/mL |
| Volume | 200 µL |
| Label | Unlabeled |
| Specificity | The high molecular weight form of Urokinase purified from human urine. |
| Reactivity | Human. Other species are not tested. |
| Storage | Maintain at 2-8°C in undiluted aliquots for up to 12 months. Do not freeze. |

Application notes

Recommended dilutions

- Western Blotting : 10 µg/mL

Other applications have not been tested.

Optimal dilutions/concentrations should be determined by the end user.

References

- 1) Husain, S. S. et al., Arch. Biochem. Biophys., 220, 31-38, 1983.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com