



Anti Prolyl-4-Hydroxylase [Clone : 6-9H6]

Product type	Monoclonal antibody
Immunogen	Rat Prolyl-4-Hydroxylase
Host	Mouse
Clone number	6-9H6
Isotype	IgG1/k
Purification	Purified
Buffer	0.1 M Na-Phosphate buffer, pH7.0 containing 0.5% protease-free bovine serum albumin and 0.05% sodium azide.
Concentration	500 µg/mL
Volume	200 µL
Label	Unlabeled
Specificity	The beta-subunit of rat Prolyl 4-Hydroxylase from the skin of newborn rat.
Reactivity	Rat. No cross-reactivity with Human. Other species are not tested.
Storage	Maintain at 2-8°C in undiluted aliquots for up to 12 months. Do not freeze.

Application notes

Recommended dilutions

- Western Blotting : 2.5 µg/mL¹⁾
- Paraffin Sections : 5 - 10 µg/mL, Trypsin treatment
- Frozen Sections : 5 - 20 µg/mL²⁾

Other applications have not been tested.

Optimal dilutions/concentrations should be determined by the end user.

References

- 1) Bai, Y. et al., J. Biochem., 99, 1563-1570, 1986.
- 2) Okuno, M. et al., J. Gastroenterol. Hepatol., 6, 271-277, 1991.



COSMO BIO CO., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com