



Anti Laminin [Clone : HL-4H3]

Product type	Monoclonal antibody
Immunogen	Human Laminin
Host	Mouse
Clone number	HL-4H3
Isotype	IgG1/k
Purification	Purified
Buffer	0.1 M Na-Phosphate buffer, pH7.0 containing 0.5% protease-free bovine serum albumin and 0.05% sodium azide.
Concentration	500 µg/mL
Volume	200 µL
Label	Unlabeled
Specificity	The P-1 fragment of human Laminin.
Reactivity	Human. Other species are not tested.
Storage	Maintain at 2-8°C in undiluted aliquots for up to 12 months. Do not freeze.

Application notes

Recommended dilutions

- Western Blotting : 5 µg/mL²⁾
- Paraffin Sections : 5 - 10 µg/mL, Trypsin treatment
- Frozen Sections : 5 - 10 µg/mL
- EIA : available²⁾

Other applications have not been tested.

Optimal dilutions/concentrations should be determined by the end user.

References

- 1) Fujimoto, N. et al., J. Med. Pharm. Sci., 23, 1273-1278, 1990 (In Japanese).
- 2) Fujimoto, N. et al., Clin. Chim. Acta, 191, 211-220, 1990.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com