



**Recombinant protein**

*For research use only, Not for diagnostic use*

**Catalog No. KN-TOYU-M04**

# Amyloid beta peptide 42 (A $\beta$ 42)

<b>Product type</b>	Recombinant Protein / Amyloid beta peptide 42 (A $\beta$ 42)
<b>Sequence</b>	DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVVIA
<b>Source</b>	E.coli
<b>Form</b>	Lyophilized
<b>Volume</b>	1 mg / tube
<b>Storage Conditions</b>	Store at 4°C. After reconstituted, store no longer than 1-2 days at 4°C, unless preservative is added. For long term storage, aliquot and store reconstituted product frozen at -20°C or lower. Aliquot to avoid cycles of freeze/thaw.

*For research use only, Not for diagnostic use.*



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9613

FAX : +81-3-5632-9619