

# Recombinant bTIMP-1

Cat. No. Y-10

Updated on Jan. 5th, 2024

## 【 I 】 Product Information

This is a purified recombinant bovine TIMP-1 from CHO cells transfected with bTIMP-1 cDNA.

Storage : -20°C

| Description         | Quantity | Concentration | Label     | Form                                                                 |
|---------------------|----------|---------------|-----------|----------------------------------------------------------------------|
| Recombinant bTIMP-1 | 100 µg   | 200 µg/mL     | Unlabeled | PBS(-): 10 mM Sodium phosphate buffer, pH7.3 containing 0.15 M NaCl. |

## 【 II 】 References

[1] Freudenstein, J. et al., mRNA of bovine tissue inhibitor of metalloproteinase: sequence and expression in bovine ovarian tissue. *Biochem. Biophys. Res. Commun.*, 171, 250-256, 1990.

**【III】 Related products**

| Description             | Cat. No. | Quantity |
|-------------------------|----------|----------|
| Human Type I Collagen   | Y-1      | 1 MG     |
| Human Type III Collagen | Y-2      | 1 MG     |
| Human Type IV Collagen  | Y-3      | 150 µG   |
| Human Type V Collagen   | Y-4      | 150 µG   |
| Human Type VI Collagen  | Y-5      | 10 µG    |
| Recombinant hTIMP-1     | Y-8      | 100 µG   |
| Recombinant hTIMP-2     | Y-9      | 100 µG   |
| Recombinant bTIMP-1     | Y-10     | 100 µG   |
| Recombinant hTIMP-1     | Y-11     | 10 µG    |
| Recombinant hTIMP-2     | Y-12     | 10 µG    |
| Recombinant bTIMP-1     | Y-13     | 10 µG    |



【JAPAN】  
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>



【Outside Japan】  
2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobiousa.com](mailto:info@cosmobiousa.com)  
URL: [www.cosmobiousa.com](http://www.cosmobiousa.com)  
Phone/FAX: (+1) 760-431-4600