

Recombinant hTIMP-1

Cat. No. Y-11

Updated on Jan. 5th, 2024

【 I 】 Product Information

This is a purified recombinant human TIMP-1 from CHO cells transfected with hTIMP-1 cDNA.

Storage : -20°C

Description	Quantity	Concentration	Label	Form
Recombinant hTIMP-1	10 µg	200 µg/mL	Unlabeled	PBS(-): 10 mM Sodium phosphate buffer, pH7.3 containing 0.15 M NaCl.

【 II 】 References

- [1] Docherty, A. J. P. et al., Sequence of human tissue inhibitor of metalloproteinases and its identity to erythroid-potentiating activity. *Nature*, 318, 66-69, 1985.
- [2] Hayakawa, T. et al., Cell growth-promoting activity of tissue inhibitor of metalloproteinases-2 (TIMP-2). *J. Cell Sci.*, 107, 2373-2379, 1994.
- [3] Kumagai, K. et al., Inhibition of matrix metalloproteinases prevents allergen-induced airway inflammation in a murine model of asthma. *J. Immunol.*, 162, 4212-4219, 1999

【 III 】 Related products

Description	Cat. No.	Quantity
Human Type I Collagen	Y-1	1 MG
Human Type III Collagen	Y-2	1 MG
Human Type IV Collagen	Y-3	150 µG
Human Type V Collagen	Y-4	150 µG
Human Type VI Collagen	Y-5	10 µG
Recombinant hTIMP-1	Y-8	100 µG
Recombinant hTIMP-2	Y-9	100 µG
Recombinant bTIMP-1	Y-10	100 µG
Recombinant hTIMP-1	Y-11	10 µG
Recombinant hTIMP-2	Y-12	10 µG
Recombinant bTIMP-1	Y-13	10 µG



【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
URL: www.cosmobiousa.com
Phone/FAX: (+1) 760-431-4600