

Human Type III Collagen

Cat. No. Y-2

Updated on Jan. 5th, 2024

【 I 】 Product Information

This is a purified human Type III Collagen.

Storage : -20°C

Description	Quantity	Concentration	Label	Form
Type III Collagen (Human)	1 mg	1 mg/mL	Unlabeled	5 mM Acetic acid

【 II 】 References

[1] Rojkind, M. et al., Increase in type I and type III collagens in human alcoholic liver cirrhosis. Proc. Natl. Acad. Sci. (USA), 73, 539-543, 1976.

【 III 】 Related products

Description	Cat. No.	Quantity
Human Type I Collagen	Y-1	1 MG
Human Type III Collagen	Y-2	1 MG
Human Type IV Collagen	Y-3	150 µG
Human Type V Collagen	Y-4	150 µG
Human Type VI Collagen	Y-5	10 µG
Recombinant hTIMP-1	Y-8	100 µG
Recombinant hTIMP-2	Y-9	100 µG
Recombinant bTIMP-1	Y-10	100 µG
Recombinant hTIMP-1	Y-11	10 µG
Recombinant hTIMP-2	Y-12	10 µG
Recombinant bTIMP-1	Y-13	10 µG

**COSMO BIO Co., LTD.****【JAPAN】**

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>

**COSMO BIO USA****【Outside Japan】**

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
URL: www.cosmobiousa.com
Phone/FAX: (+1) 760-431-4600