

## 抗ヒトAFPモノクローナル抗体

品番	KR-006
クローン	1D5
型	IgG1, $\kappa$
由来	BALB/Cマウス腹水
精製法	硫酸塩析, r-プロテインA
溶液	PBS (0.09%アジ化ナトリウム含む)
数量	100 $\mu$ g
保存法	4~8°C 凍結融解は出来るだけ避けて下さい。
備考	抗ヒトAFP抗体(1D5)はサンドイッチELISA法の二次抗体に適しています。なお、固相としては抗ヒトAFP抗体(6D2)が適しています。



人と科学のステキな未来へ  
**コスモ・バイオ株式会社**

— 商品に関するお問い合わせ —

TEL: 03-5632-9610 FAX: 03-5632-9619

本社所在地 〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル

**MOUSE ANTI AFP  
MONOCLONAL ANTIBODY**

**Catalog No.** LNM-KR-006

**CLONE:** 1D5

**QUANTITY:** 0.1mg

**SPECIFICITY:** Human  $\alpha$  fetoprotein

**IMMUNOGEN:** Human  $\alpha$  fetoprotein

**Type/Subtype:** IgG1  $\kappa$

**APPLICATIONS:** Sandwich ELISA (This antibody can be used as a detection antibody in a human  $\alpha$  fetoprotein ELISA combination with the monoclonal capture antibody 6D2.)

**FORMAT:** Purified immunoglobulin from BALB/C Mouse Ascites  
(Affinity chromatography)

**PRESENTATION:** Liquid in 15mM Phosphate buffer containing 0.15M NaCl,  
0.09% sodium azide, pH 7.2

**STORAGE/HANDLING:** 2 to 8°C, Avoid repeated freeze/thaw cycle.

***For research use only; not for use as a diagnostic.***



**COSMO BIO CO., LTD.**

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN  
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com