



## TNP-Ascaris

<b>Form</b>	Lyophilized protein powder 2,4,6-trinitrophenyl Ascaris Crude Extract
<b>Volume</b>	Protein: 5mg Conjugation ratio: 33.7 TNP groups/molecule (Lot #748021) (As average of molecular weight: 100,000)
<b>Concentration</b>	Add 0.5mL of distilled water. It will become 5mM Borate buffered saline [BBS] (pH 8.0), 0.9% NaCl solution.
<b>Storage</b>	Store below 4°C (below -20°C for prolonged storage or after dissolved)

### Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG
LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

### Adjuvant

LG-6000	ALUM
---------	------

*For research use only. Not for clinical diagnosis.*



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

## TNP-Ascaris

## I. 内容

lot No. 748021


性 状	2,4,6-trinitrophenyl Ascaris Extract 凍結乾燥粉末
規 格	蛋白質量：5mg TNP基数：1分子あたり33.7個（平均分子量 100,000としての計算値） 0.5mlの精製水を加えることで5mM硼酸緩衝溶液 （pH8.0）0.9%NaCl溶液になります。

## II. 保存上の注意

4℃以下（長期間又は溶液の場合は-20℃以下）で保存して下さい。

製造元

総発売元

 株式会社 エル・エス・エル

コスモ・バイオ株式会社

## 新製品案内

LG-6061	OA Injection FFDP
LG-6202	EB Injection FFDP
LG-5779	Hinoki Pollen Extract
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LB-8109	Rat Anti TNP-Ascaris