



For research use only. Not for clinical diagnosis.

Catalog No. 2ZTHBR2

THBR2

Code	2ZTHBR2
Description	Mouse IgM, Purified
Source	Mouse ascites
Application	<p>This product inhibits nonspecific reaction caused by human anti-mouse antibody (HAMA) in immunoassays using mouse monoclonal antibodies.</p> <p>The recommended concentration is 0.01-0.1 mg/mL for ELISA, however, the optimum condition should be determined for each assay.</p>
Purification	Ammonium sulfate precipitation
Subclass	IgM
Form	Saline (containing 0.1% NaN ₃)
Protein Concentration	10mg/mL
Volume	2mL
Storage	<p>Maintain at 2-8 °C (Below -20 °C for long term storage)</p> <p>Avoid repeated freeze/thaw cycles</p>



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619



人と科学のステキな未来へ

コスモ・バイオ株式会社

For research use only. Not for clinical diagnosis.

Catalog No. 2ZTHBR2

THBR2

Code	2ZTHBR2
Description	マウス IgM 精製品
Source	マウス腹水
Application	マウス抗体を用いたイムノアッセイにおける、マウスに対する異好性抗体 (HAMA:Human Anti-Mouse Antibody) による非特異反応を抑制します。 サンドイッチ測定系における標準添加濃度：10～100 μ g/mL あらかじめ各測定系での最適条件を検討してください。
Purification	硫酸アンモニウム塩析精製
Subclass	IgM
Form	0.1%アジ化ナトリウムを含む生理的食塩水溶液
Protein Concentration	10mg/mL
Volume	2mL
Storage	2～8°C 3ヶ月 -20°C 3年間 (凍結融解の繰り返しは避けてください)



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619