





GCGR Monoclonal Antibody

Product Code	CSB-RA009316A1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47871
Immunogen	Recombinant Human GCGR protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer; Metabolism; Signal transduction
Gene Names	GCGR
Accession NO.	9G1

Image



The Binding Activity of Human GCGR with Anti-GCGR recombinant antibody.
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human GCGR (CSB-MP009316HU1) at 2 µg/mL can bind Anti-GCGR recombinant antibody, the EC₅₀ is 14.57-32.56 ng/mL.



The Binding Activity of Human GCGR with Anti-GCGR recombinant antibody.
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human GCGR (CSB-YP009316HU1) at 2 µg/mL can bind Anti-GCGR recombinant antibody, the EC₅₀ is 104.1-179.5 ng/mL.