



# S Antibody

<b>Product Code</b>	CSB-RA33245A0GMY
<b>Abbreviation</b>	S
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0DTC2
<b>Immunogen</b>	Recombinant Human Novel Coronavirus Spike glycoprotein (S) (16-685aa)
<b>Species Reactivity</b>	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Alias</b>	S, S1, S1-RBD, Spike glycoprotein
<b>Species</b>	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
<b>Research Area</b>	Microbiology
<b>Gene Names</b>	S