



# SARS-CoV-2 Spike RBD Nanobody, HRP conjugated

<b>Product Code</b>	CSB-RA33245B2GMY
<b>Abbreviation</b>	S
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0DTC2
<b>Immunogen</b>	Recombinant Human Novel Coronavirus Spike glycoprotein(S) (319-541aa) (CSB-YP3324GMY1 and CSB-MP3324GMY1b1)
<b>Species Reactivity</b>	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
<b>Tested Applications</b>	ELISA; Recommended dilution: ELISA:1:1000-1:500000
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	VHH fusion with human IgG1 Fc
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Recombinant Antibody
<b>Immunogen Species</b>	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
<b>Research Area</b>	Microbiology
<b>Accession NO.</b>	A1