



N Antibody

Product Code	CSB-RA33255A0GMY
Abbreviation	N
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DTC9
Immunogen	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa)
Species Reactivity	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG1
Clonality	Monoclonal
Alias	Nucleocapsid protein, NC, protein N, N
Species	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
Research Area	Microbiology
Gene Names	N