

PEPTIDE

His tag peptide

Code No.	Quantity	Form
3310-205	2 mg x 5 vials	Lyophilized

DESCRIPTION:

This synthetic peptide including 6xHis tag sequence has been designed for the elution of His tagged Protein PURIFICATION KIT. The peptide sequence is XXX-(6xHis)-XXX.

MOLECULAR WEIGHT:

1,616 Da

SUPPLIED:

Lyophilized in PBS.

Reconstitute in 1.0 mL of distilled water. After reconstitution, the His tag peptide concentration will be 2 mg/mL in PBS.

STORAGE:

Lyophilized His tag peptide is stable for one year from the date of shipment when stored at 4°C. After reconstitution, prepare appropriate aliquots and freeze them at -20°C. Avoid repeated freeze/thaw cycles.

RELATED PRODUCTS:

Please visit our website at <https://ruo.mbl.co.jp/>.

220628-4