

Anti-Mincle (Mouse) mAb

(Functional Grade)

CODE No.	D292-3M2
CLONALITY	Monoclonal
CLONE	4A9
ISOTYPE	Rat IgG1 κ
QUANTITY	100 μ L, 1 mg/mL
SOURCE	Purified IgG from hybridoma supernatant
IMMUNOGEN	RBL-2H3 cells expressing full-length mouse Mincle
REACTIVITY	This clone reacts with mouse Mincle (Clec4e).
FORMULATION	PBS. Azide free, 0.22 μ m sterile-filtered Endotoxin level is < 0.5 EU/mg antibody, as determined by the LAL assay.
STORAGE	This antibody solution is stable for one year from the date of purchase when stored at -20°C.
APPLICATION	
<u>Functional activity</u>	1-10 μ g/mL for blocking, Reference 1)
Entrez Gene ID	56619 (Mouse)
REFERENCES	1) Yamasaki, S., <i>et al.</i> , <i>Nat. Immunol.</i> 9 , 1179-1188 (2008) 2) Yamasaki, S., <i>et al.</i> , <i>PNAS</i> 106 , 1897-1902 (2009) 3) Ishikawa, E., <i>et al.</i> , <i>J. Exp. Med.</i> 206 , 2879-2888 (2009)

For more information, please visit our web site <https://ruo.mbl.co.jp/>.

200310-2